



**PORTABLE SELF-ASSESSMENT AUDIOMETER
USING RASPBERRY PI**

MINI PROJECT REPORT

Submitted by

NAME	ROLL NO
ADHIL.M	727622BEC114
PRANESH .S	727622BEC018
NAVEEN.S	727622BEC110

Under the guidance of

Mr.A.Shafeek,M.E.,

In partial fulfillment off the requirement for the award of

BACHELOR OF ENGINEERING

In

ELECTRONICS AND COMMUNICATION

ENGINEERING

Dr. MAHALINGAM COLLEGE OF

ENGINEERING AND

TECHNOLOGY

An Autonomous Institution

Affiliated to ANNA UNIVERSITY

CHENNAI – 600 025

MAY-2024

**Dr. MAHALINGAM COLLEGE OF ENGINEERING AND
TECHNOLOGY, POLLACHI -642 003**

(An Autonomous Institution Affiliated to Anna University, Chennai - 600 025)

BONAFIDE CERTIFICATE

Certified that this Mini project report titled "**PORTABLE SELF-ASSESSMENT AUDIOMETER USING RASPBERRY PI**" is the bonafide work of

NAME	ROLL NO
ADHIL.M	727622BEC114
PRANESH .S	727622BEC018
NAVEEN.S	727622BEC110

who carried out the mini project under my supervision.

Dr.R.Sudhakar,M.E.,Ph.D

Mr.A.Shafeek M.E., AP(SS)

HEAD OF THE DEPARTMENT

Department of ECE
Dr. Mahalingam College of
Engineering and Technology
Pollachi- 642003

SUPERVISOR

Department of ECE
Dr. Mahalingam College of
Engineering and Technology
Pollachi- 642003

Submitted for the Autonomous End Semester Mini Project Examination

held on _____

INTERNAL EXAMINER

EXTERNAL EXAMINER

**Dr. MAHALINGAM COLLEGE OF ENGINEERING AND
TECHNOLOGY, POLLACHI -642 003**

(An Autonomous Institution Affiliated to Anna University, Chennai - 600 025)

TECHNOLOGY READINESS LEVEL (TRL) CERTIFICATE

Project Title: Portable Self-Assessment Audiometer Using Raspberry Pi

Course Code: 19ECPN6401- Mini Project

NAME	ROLL NO
ADHIL.M	727622BEC114
PRANESH .S	727622BEC018
NAVEEN.S	727622BEC110

Guide Name: Mr.A.Shafeek M.E., AP(SS)/ECE

Technology Readiness Level* (TRL) of this Project: TRL 06

Prototype System Verified (System/process prototype demonstration in an operational environment)

SIGNATURE OF THE GUIDE

HOD

INTERNAL EXAMINER

EXTERNAL EXAMINER

(*Refer to TRL Definition and for projects, Minimum TRL should be 4)

ACKNOWLEDGEMENT

We wish to express our sincere thanks to all who have contributed to do this project through their support, encouragement and guidance.

We extend our gratitude to our management for having provided us with all facilities to build my project successfully. We express our sincere thanks to our honourable Secretary, **Dr.C.Ramaswamy, M.E., Ph.D., F.I.V.**, for providing the required amenities.

We take this opportunity to express our deepest gratitude to our principal **Dr.P.Govindhasamy, M.Tech., Ph.D.**, who provide suitable environment to carry out the project.

We express our extreme gratefulness to **Dr.R.Sudhakar, M.E., Ph.D.**, Professor and Head of the department of Electronics and Communication Engineering, for his constant support and encouragement.

We wish to express our deep sense of gratitude and thankfulness to my project guide **Mr.A.Shafeek,M.E.**, for the valuable suggestion and guidance offered during the course of the project.

Finally, we are committed to place our heartfelt thanks to all those who had contributed directly and indirectly towards the success of the completion of this project.

PORTABLE SELF-ASSESSMENT AUDIOMETER

USING RASPBERRY PI

ABSTRACT

People with hearing impairments often face irreversible damage. To determine the extent of hearing loss across different frequencies in each ear, a hearing screening test is used. However, traditional audiometers require an audiologist to conduct the test, which can be time-consuming and expensive for the individual. This study aims to create a new portable audiometer for self-assessment of hearing.

The portable audiometer consists of a computer, Raspberry Pi 3 B+, patient response button, and headphones. Sound signals are delivered to the patient through the headphones, and the patient responds using the left mouse button. Based on the patient's responses, an automatic audiogram is generated, showing the relationship between frequency and intensity, which indicates the volume of sound pressure.

The results, including the audiogram and raw data, are saved in CSV files named with the time and date of the test. A familiarization procedure is employed to help the hearing-impaired individual understand and perform the response task. The efficient Hughson Westlake procedure, which is less time-consuming, is implemented in Python, a popular open-source programming language, to obtain the audiogram. Using Python helps reduce software development costs.

BATCH 08

ADHIL.M - 727622BEC114

PRANESH .S - 727622BEC018

NAVEEN.S - 727622BEC110

SIGNATURE OF THE SUPERVISOR

Mr.A.Shafeek M.E., AP(SS)/ECE

Department of ECE

Dr. Mahalingam College of

Engineering and Technology

Pollachi- 64200

TABLE OF CONTENTS

S.NO	CHAPTER	PG.NO
	INTRODUCTION	I
	CERTIFICATE	II
	TECHNOLOGY READINESS LEVEL (TRL) CERTIFICATE	III
	ACKNOWLEDGEMENT	IV
	ABSTRACT	V
	TABLE OF CONTENT	VI
	LIST OF FIGURES	IX
	LIST OF TABLES	XII
	LIST OF ABBRIVATIONS	XIII
1	INTRODUCTION	1
	1.1 INTRODUCTION	1
	1.2. MOTIVATION	3
	1.3 OVERVIEW	3
	1.4 OBJECTIVES	3
	1.5 PROBLEM STATEMENT	4
	1.6. BENEFITS OF THE PROJECT	4
2	LITERATURE REVIEW	5
	2.1 INTRODUCTON	5
3	HEARING LOSS	8
	3.1 INTRODUCTION	8
	3.2 ANATOMY OF THE EAR	8
	3.3. PREVALENCE OF HEARING LOSS	9
	3.4 CAUSES AND RISK FACTORS	9
	3.5. TYPES OF HEARING LOSS	10
	3.6. SYMPTOMS OF HEARING LOSS	11
	3.7 STAGES OF HEARING LOSS	11

4	SELF-ASSESSMENT AUDIOMETER	12
4.1	INTRODUCTION	13
4.2.	PORTABLE AUDIOMETER UTILIZING RASPBERRY PI	13
4.3.	PORTABLE AUDIOMETER (RASPBERRY PI)	14
	VS MOBILE HEARING TEST	
4.4.	SELF-ASSESSMENT AUDIOMETER FOR TESTING:	15
	POSSIBLE AND NOT POSSIBLE SCENARIOS	
4.5.	EXISTING VS PROPOSED SELF-ASSESSMENT AUDIOMETERS	16
5.	PURE TONE AUDIOMETER	18
5.1.	INTRODUCTION	18
5.2.	PURPOSE OF PTA	18
5.3.	TYPES OF TESTING	19
5.4.	AUDIOGRAM	20
5.5.	AUDIOGRAM CONFIGURATION	22
5.6.	LIMITATIONS	22
5.7.	APPLICATION	23
5.8.	FEATURES	23
6	IMPLEMENTATION OF HARDWARE MODULES	24
6.1	INTRODUCTION	24
6.2.	BLOCK DIAGRAM	24
	COMPONENT DESCRIPTION	
6.3.	RASPBERRY PI	25
	6.2.1. FEATURES OF RASPBERRY PI	25
	6.2.2. RASPBERRY PI GPIO'S PIN	26
	6.2.3. RASPBERRPY PI PIN DETAILS	27
6.4.	RASPBERRY PI CASE	28
6.5.	MOUSE	28
6.6.	HEADPHONE	29
	6.6.1. TDH 49 HEADPHONE CHARACTERISTICS AND USAGE	30
6.7.	ETHERNET CABLE	30
6.8.	POWER SUPPLY	31

7	SOFTWARE IMPLEMENTATION	32
	7.1 INTRODUCTION	32
	7.2. RASPBERRY PI OS (RASPBIAN)	32
	7.3. VNC: REMOTE ACCESS A RASPBERRY PI	33
	7.4. PYTHON ON RASPBERRY PI: THONNY EDITOR	34
8.	VALIDATION AND FEEDBACK: MEETING WITH ENT DOCTOR	35
	8.1. INTRODUCTION	35
	8.2. DISCUSSION WITH DR. B. HARSHA	35
	8.3. RECOMMENDATION TO MEET AN AUDIOLOGIST	36
	8.4. DISCUSSION WITH THE ENT DOCTOR	36
	8.5. TESTING THE SELF-ASSEMBLED PURE TONE AUDIOMETER	37
	8.6. FEEDBACK AND CONCLUSION	37
9.	MEETING WITH AUDIOLOGIST BALACHANDRAN	38
	9.1. INTRODUCTION	38
	9.2. MEETING WITH AUDIOLOGIST BALACHANDRAN	38
	9.3. PROJECT DISCUSSION AND RECOMMENDATIONS	39
	9.4. QUESTIONS AND AUDIOMETER INSIGHTS	39
	9.5. GUIDANCE AND LEARNING	40
	9.6. KEY FEEDBACK FROM AUDIOLOGIST	40
	9.7. EXPLORATION AND FUTURE SCOPE	41
	9.8. CONCLUSION	41
10	RESULTS	42
	10.1 INTRODUCTION	42
	10.2 TESTING OF FINAL WORKING MODEL	42
	10.3 WORKING MODEL	43
	10.4 RESULT ANALYSIS	44
11	CONCLUSION	48
12	BIBLIOGRAPHY	49
13	APPENDICES	51
	PLAGIARISM REPORT	58

LIST OF FIGURES

S NO	FIGURES	TOPIC	PAGE NO
1.	1.1	PURE TONE AUDIOMETRY (PTA)	2
2.	3.2	ANATOMY OF THE EAR	8
3.	3.5	TYPES OF HEARING LOSS	10
4.	3.6	SYMPTOMS OF HEARING LOSS	11
5.	3.7	STAGES OF HEARING LOSS	11
6.	4.1	SELF-ASSESSMENT AUDIOMETER	13
7.	4.3.1	PORTABLE AUDIOMETER	14
8.	4.3.2	MOBILE HEARING TEST	14
9.	4.5.1	PORTABLE SELF ASSESSMENT AUDIOMETER	16
10..	4.5.2.	AUDIOMETER	16
11.	5.1.	PURE TONE AUDIOMETER	18
12	5.3.	TYPES OF TESTING	19
13	5.4.	AUDIOGRAM	20
14.	5.4.1	Audiogram Sheet	21

15.	5.5.	AUDIOGRAM CONFIGURATION	22
16.	6.2.	Block Diagram	24
17	6.3.1	RASPBERRY PI	25
18.	6.3.2	PIN DESCRIPTION OF RASPBERRY PI	26
19.	6.3.3	RASPBERRY PI GPIO'S PIN LAYOUT	27
20.	6.3.4	RASPBERRY PI GPIO	27
21.	6.4	RASPBERRY PI CASE	28
22.	6.5	MOUSE	29
23.	6.6	HEADPHONE	29
24.	6.7	ETHERNET CABLE:	31
25.	6.8	POWER ADAPTER	31
26.	7.2	INSTALLING RASPBIAN.	32
27.	7.2.1	THE RASPBERRY PI DESKTOP.	33
28.	7.3	VNC: REMOTE ACCESS TO A RASPBERRY	34
29.	7.4	THONNY IDE	34
30.	8.2.	CONSULTING WITH THE DOCTOR	41

31.	9.2.	CONSULTING WITH THE AUDIOLOGIST	44
32	9.3.	AUDIOMETRY TEST	45
33.	9.6.	AUDITIVO AUDIOMETER	46
34.	10.3.	OVERALL INITIALIZATION PURE TONE AUDIOMETER	50
35.	10.4.1.	PYTHON CODE IN THONNY IDE	50
36.	10.4.2.	INSTRUCTION OF THE PROJECT	51
37.	10.4.3.	PATIENT RESPONSE OF THE AUDIOMETER	51
38.	10.4.4.	OUTPUT OF THE AUDIOGRAM CHART	52
39	10.4.5.	RESULT AUDIOGRAM CHART	52
40	10.4.6.	PORTABLE SELF ASSESSMENT AUDIOMETER	53
41	11	PLAGIARISM REPORT	55

LIST OF TABLES

S NO	TABLE	TOPIC	PAGE NO
1.	3.7	Stages of hearing loss	12
2.	4.3	Portable Audiometer (Raspberry Pi) Vs Mobile Hearing Test	15
3.	4.7	Existing vs proposed self-assessment audiometers	16
4,	10.4	Result analysis	53

LIST OF ABBRIVATIONS

1. HA - Hearing aid
2. SNHL - Sensorineural hearing loss
3. WNL - Within normal limits
4. AU - Both sides (ears)
5. AS - Left
6. AD - Right
7. BC - Bone conduction
8. AC - Air conduction
9. PTA - Pure-tone average
10. UCL - Uncomfortable loudness level
11. MCL - Most comfortable loudness level
12. HFA - High frequency average
13. HL - Hearing level
14. SRT - Speech reception threshold
15. SAT - Speech awareness threshold

CHAPTER - 1

INTRODUCTION

1.1 INTRODUCTION

Hearing impairment is a widespread chronic condition affecting adults globally, especially as they age, particularly at frequencies exceeding 2000 Hz. According to the World Health Organization (WHO), approximately 466 million people worldwide are affected by hearing loss, with 75% residing in developing nations. Factors such as prolonged exposure to loud noises, aging, tumors, illnesses, and ototoxic medications contribute to hearing loss, which may occur alone or alongside tinnitus. Though not fatal, hearing loss can lead to depression, communication difficulties, decreased functional abilities, and social isolation. Detecting and treating hearing loss early is crucial to enhance the quality of life for older adults dealing with hearing impairment.

Research conducted in the United Kingdom found that individuals who identified their hearing loss early and consistently used hearing aids experienced greater benefits compared to those who delayed using hearing aids until later stages of the condition. Therefore, it is essential to conduct hearing screening tests in the early stages to minimize the impact of hearing loss on individuals' lives. Hearing levels are typically evaluated across frequency ranges centered at 125Hz, 250 Hz, 500 Hz, 1 kHz, 2 kHz, 4 kHz, and 8 kHz. Although the human hearing range is commonly stated as 20 to 20,000 Hz, variations can occur, especially at higher frequencies.

Audiometers are instruments used to assess hearing capacity, often employing pure tone audiometry (PTA) testing, where individuals respond to presented pure tones. Traditionally, audiologists administer these tests in soundproof booths, necessitating patient visits to clinics, which may be inconvenient, especially for elderly individuals with mobility issues. During the test, patients wear sound-isolating headphones, and tones ranging from 125 Hz to 8 kHz are played at different intensities to determine hearing thresholds. The intensity level is adjusted to reach the patient's threshold, following the Hughson Westlake ascending method. While traditional audiometers are precise, they are costly and time-consuming.

Advances in information and communication technology have led to the development of screening audiometers and portable hearing-test devices, offering improved efficiency and

cost-effectiveness. Various approaches, such as auto kits utilizing personal computers or portable audiometers with dedicated hardware units, have been devised to streamline the testing process. However, challenges such as equipment costs or the necessity for audiologist involvement persist.

The present study introduces a self-assessment hearing test aimed at being cost-effective, straightforward, and user-friendly. The goal was to develop a computer-based screening kit that individuals with minimal computer proficiency could easily utilize to conduct their own hearing tests. This proposed pure tone audiometry (PTA) system offers high portability and affordability, with the aim of broadening access to hearing screening.



Fig 1.1 Pure Tone Audiometry (PTA)

1.2 MOTIVATION

The motivation for developing the portable self-assessment audiometer arises from the need to overcome limitations of traditional testing methods. This includes improving accessibility, reducing costs, and enhancing convenience, especially for those in remote areas. By leveraging technology and affordability, the project aims to democratize hearing healthcare, empowering individuals to monitor their hearing health easily and effectively. The device's user-friendly design and data storage capabilities further support its role in facilitating early detection and intervention, ultimately improving outcomes and quality of life for individuals with hearing impairments.

1.3 OVERVIEW

The project focuses on developing a Portable Audiometer for Hearing Self-Assessment, aiming to address the challenges posed by traditional audiometry testing in terms of accessibility, cost, and time consumption. Leveraging Raspberry Pi, Python programming language, and audio hardware, the audiometer enables individuals to assess their hearing independently, offering a cost-effective and user-friendly solution.

The audiometer comprises hardware components including Raspberry Pi 3 B+, a mouse for user interaction, and headphones for audio output. Raspberry Pi serves as the core computing platform, generating pure tone signals across different frequencies and intensity levels. Python programming language is utilized for software development, managing signal generation, user interaction, response recording, and data analysis.

During the testing procedure, the audiometer conducts a series of tests for both ears, playing pure tone signals at varying frequencies and intensity levels. Users respond to the audible tones by clicking the mouse button. The modified Hughson Westlake procedure is employed to adjust signal intensity based on user response

1.4 OBJECTIVES

The project endeavors to create a Portable Audiometer for Hearing Self-Assessment, leveraging Raspberry Pi, Python programming, and audio hardware. Its goal is to address the drawbacks of traditional audiometry by providing a budget-friendly, easy-to-use solution for self-directed hearing evaluation. By producing pure tone signals and integrating a modified version of the Hughson Westlake procedure, individuals can react to audible tones by clicking the mouse.

The gathered data will undergo analysis to ascertain hearing thresholds and produce audiograms, visually depicting hearing levels. Ultimately, the aim is to improve the accessibility of hearing screening, facilitating early identification and management of hearing impairment

1.5. PROBLEM STATEMENT

Traditional audiometric testing methods present barriers to individuals seeking to assess their hearing health, including limited accessibility, high costs, and inconvenience. Particularly challenging for those in remote areas, these methods often result in delays in diagnosis and intervention. To address these issues, there is a pressing need for a portable self-assessment audiometer. This device must be affordable, user-friendly, and capable of providing accurate assessments without professional oversight. Additionally, it should offer convenient data storage and analysis for collaboration with healthcare professionals. The challenge is to design and implement a portable audiometer that democratizes access to hearing healthcare, empowering individuals to monitor and manage their hearing health effectively. It must be cost-effective, portable, easy to use, and capable of delivering accurate results, ultimately improving outcomes for individuals with hearing impairments.

1.6. BENEFITS OF THE PROJECT

The project of developing a portable self-assessment audiometer offers several benefits:

1. **Accessibility:** It provides easy access to hearing assessments for individuals, especially in remote areas where audiologists are scarce.
2. **Cost-Effectiveness:** By utilizing affordable components and open-source software, it reduces the financial burden associated with traditional audiometric tests.
3. **Convenience:** Users can conduct tests at their convenience, avoiding the need for scheduling appointments and traveling to clinics.
4. **Early Detection:** Regular self-assessment allows for the early detection of hearing issues, leading to timely intervention and improved outcomes.
5. **Empowerment:** It empowers individuals to take control of their hearing health by facilitating regular monitoring and self-management.
6. **Portability:** The device's portability enables testing in various settings, including homes, schools, workplaces, and community centers.

CHAPTER -2

LITERATURE REVIEW

2.1 INTRODUCTION

This research paper addresses the need for accessible hearing health assessment methods. Traditional audiometry faces barriers like cost and accessibility. To overcome these, we explore a portable audiometer for self-assessment, utilizing Raspberry Pi technology. This innovation aims to empower individuals to monitor their hearing health independently, particularly benefiting underserved communities. Through rigorous testing, we aim to demonstrate the effectiveness of this solution in improving accessibility and inclusivity in hearing healthcare.

2.2. METHODOLOGY

[1] “A NOVEL RASPBERRY PI-3 BASED PURE TONE AUDIOMETER AND VERIFICATION OF CALIBRATION WITH STANDARD SYSTEM”

M.Dharani Kumar Chowdhary , Dr. C. Nagaraja , Dr. C. Sandeep Kumar Reddy , Dr.Srinivasarao Udara ,worked on “A novel raspberry pi-3 based pure tone audiometer and verification of calibration with standard system” In the present era, The Raspberry Pi-based audiometer utilizes a sound card and headphones to generate and present pure tones at various frequencies and intensities, mimicking a traditional audiometer. Users would indicate hearing the tones through a method not specified in the current summary (possibly a button press). A crucial aspect of this system is calibration, ensuring accurate results.

The study emphasizes verifying the Raspberry Pi system's calibration against a standard audiometer. This portable and potentially lower-cost solution holds promise for increasing access to hearing screenings in areas with limited resources. It might also be useful for individuals to monitor their hearing health at home, although a professional diagnosis remains essential.

[2] "DEVELOPMENT OF HEARING SELF-ASSESSMENT PURE TONE AUDIOMETER"

Marwa Gargouri, Mondher Chaoui, and Patrice Wira present a pioneering study on the "Development of hearing self-assessment pure tone audiometer." Addressing the global prevalence of hearing loss and its consequential impact, the researchers propose an innovative solution to provide accessible and cost-effective hearing screening. Their portable audiometer integrates Raspberry Pi technology with a computer, patient response button, and headphones, allowing individuals to conduct self-administered hearing tests. By employing the Hughson Westlake procedure implemented in Python, the audiometer automatically generates audiograms based on patient responses. This approach aims to mitigate the time and cost burdens associated with traditional audiology, particularly benefiting elderly and underserved populations. The study underscores the significance of early detection and intervention in managing hearing impairment while advancing accessibility to screening tools through innovative technology and open-source software.

[3] HEARTHAT? - AN APP FOR DIAGNOSING HEARING LOSS

Silvia Figueira, Kevin Nguyen, and Shweta Panditrao's project, "HearThat? - An app for diagnosing hearing loss," showcased the transformative potential of smartphone applications in tackling healthcare obstacles. Their work aligns with the essence of our project, focusing on leveraging technology to enhance accessibility and efficiency in healthcare. By developing an app specifically designed for diagnosing hearing loss, they exemplify the innovative possibilities of mobile health solutions. This project underscores the relevance of incorporating smartphone-based tools into healthcare practices, emphasizing the importance of accessibility and ease of use in addressing medical challenges. Their insights contribute to the broader landscape of utilizing technology to empower individuals and improve healthcare outcomes.

[4] "PORTABLE AUDIOMETER FOR DETECTING HEARING DISORDER AT AN EARLY STAGE FOR CANCER PATIENTS

In their work presented at the 2016 International Conference on Automatic Control and Dynamic Optimization Techniques, Ritu Rani and H.T. Patil introduced a "Portable audiometer for detecting hearing disorder at an early stage for cancer patients." This innovative device signifies a multifaceted approach to healthcare, emphasizing the significance of early detection

in managing health conditions. By focusing on cancer patients, the research addresses the often-overlooked aspect of hearing impairment resulting from cancer treatments or related complications. The development of a portable audiometer tailored to this specific demographic underscores the importance of personalized medical solutions and the integration of technology into healthcare delivery. This project likely involved the design and validation of a specialized diagnostic tool capable of early detection, thereby highlighting the intersection of medical research, technology, and patient care.

[5] The conception of a software pure tone audiometer application by Marwa Gargouri, Mondher Chaoui, and Abdennaceur Kachouri represents a significant milestone in audiology research and technological innovation. Their work embodies a shift towards digital solutions for conducting pure tone audiometry, offering the potential to revolutionize the field by improving accessibility and reducing costs associated with traditional audiometers. By leveraging software and technology, their application not only enhances the efficiency and accuracy of hearing screening methods but also addresses the challenges faced by individuals with hearing impairments. Integrating their research findings into the literature review of your project could provide comprehensive insights into the development and application of software-based audiometry tools, further advancing the understanding and implementation of hearing assessment technologies.

[6] Fei Chen and Shuai Wang explore the development of smartphone-based self-hearing diagnosis using hearing aids. Their research investigates the potential of leveraging smartphone technology to offer accessible and convenient methods for individuals to assess their hearing abilities. By integrating hearing aid functionality with smartphone applications, the study aims to provide users with real-time feedback on their hearing health, empowering them to proactively monitor and manage any changes in their auditory capabilities. The utilization of smartphones as diagnostic tools for hearing impairment holds promise in extending the reach of hearing healthcare services to broader populations, particularly those with limited access to traditional clinical resources. Chen and Wang's work contributes to the advancement of tele audiology and mobile health solutions, fostering greater independence and autonomy in managing hearing-related issues.

CHAPTER – 3

HEARING LOSS

3.1. INTRODUCTION

Current approaches for evaluating hearing loss, like traditional audiometry, demand skilled personnel and specialized tools, rendering them slow, expensive, and out of reach for many. There's a demand for a portable and easy-to-use audiometer that empowers individuals to assess their own hearing, especially those with limited computer proficiency. Such a device should provide an affordable means for frequent monitoring of hearing health, boosting access to hearing care services and enriching the lives of those with hearing challenges.

3.2 ANATOMY OF THE EAR




Fig 3.2 Anatomy of the Ear

The ear consists of three primary components: the external ear (comprising the pinna and ear canal), middle ear (consisting of the ossicles), and inner ear (housing the cochlea and vestibular system). Sound waves are captured by the pinna, then transmitted through the ear

canal to the eardrum. The eardrum vibrates in response, conveying these vibrations through the ossicles to the cochlea. Within the cochlea, numerous hair cells transform these vibrations into electrical signals, which are then relayed to the brain via the auditory nerve, allowing for auditory perception. Additionally, the inner ear's vestibular system plays a crucial role in maintaining equilibrium. Understanding the anatomical structure of the ear is essential for diagnosing and managing conditions affecting hearing and balance.

3.3 PREVALENCE OF HEARING LOSS:

Hearing loss is widespread, with statistics indicating its significant impact across different demographics:

1. In India, over 63 million people (approximately 6% of the population) suffer from significant hearing impairment.
2. Factors contributing to this include age, gender, socioeconomic status, and environmental conditions.
3. Efforts to address this issue include awareness campaigns, early detection programs, and improvements in healthcare infrastructure.
4. There remains a significant need for continued investment in accessible and affordable hearing healthcare services.

3.4. CAUSES AND RISK FACTORS

1. Aging: Natural deterioration of auditory function over time.
2. Genetics: Inherited predispositions to hearing loss.
3. Noise Exposure: Prolonged exposure to loud noises.
4. Ototoxic Medications: Certain medications that can damage hearing.
5. Medical Conditions: Diabetes, cardiovascular disease, and ear infections.

3.5. TYPES OF HEARING LOSS




Fig 3.5. Types of Hearing Loss

1. Sensorineural Hearing Loss:

This form of hearing impairment typically arises within the cochlea or inner ear and is often irreversible. A variety of factors, such as genetic predisposition, prenatal or childhood illnesses, specific medications, accidents leading to inner ear damage, or exposure to high levels of noise, can contribute to its development.

2. Conductive Hearing loss:

Conductive hearing loss occurs when there is an obstruction or malfunction in the transmission of sound to the inner ear, while the inner ear itself remains unaffected. This type of hearing loss may be temporary or permanent. Surgical procedures or the use of hearing aids can help mitigate the blockage, allowing sound to reach the inner ear.

3. Combined Hearing loss:

Mixed hearing loss encompasses aspects of both sensorineural and conductive hearing loss. It involves challenges related to both the inner ear and the transmission of sound. Treatment approaches may involve a blend of surgical interventions, hearing aids, or other customized strategies to address both forms of impairment.

3.6. SYMPTOMS OF HEARING LOSS:




Fig 3.6. Symptoms of Hearing Loss

1. Difficulty hearing speech, especially in noise.
2. Asking for repetition frequently.
3. Increasing volume on devices.
4. Speech sounding unclear.
5. Struggling with high-pitched sounds.
6. Hearing ringing or buzzing.
7. Experiencing ear discomfort.

3.7. STAGES OF HEARING LOSS




Fig 3.7. Stages of hearing loss

- 1. Normal Auditory Function (0-25 dB):** Adequate capacity to perceive faint sounds, whispers, and gentle speech with clarity and ease.
- 2. Mild Auditory Impairment (26-40 dB):** Challenges in discerning soft or distant speech, particularly amid noisy surroundings.
- 3. Moderate Auditory Impairment (41-55 dB):** Effort required to comprehend conversational speech, especially when background noise is present.
- 4. Moderately Severe Auditory Impairment (56-70 dB):** Struggles in grasping speech without external assistance, such as hearing aids.
- 5. Severe Auditory Impairment (71-90 dB):** Limited capability to comprehend speech without external amplification; reliance on sign language may be necessary
- 6. Profound Hearing Loss (>90 dB):** Very limited or no hearing without hearing aids or cochlear implants; heavily reliant on visual cues for communication.

SL. NO.	STAGE	DESCRIPTION	THRESHOLDS (dB)
1.	NORMAL	No significant impairment	0-25
2.	MILD	Difficulty with faint or distant speech	26-40
3.	MODERATE	Trouble with conversational speech	41-55
4.	MODERATELY SEVERE	Difficulty even in quiet environments	56-70
5.	SEVERE	Significant difficulty without amplification	71-90
6.	PROFOUND	Very limited or no hearing without assistance	>90

TABLE 1: STAGES OF HEARING LOSS

CHAPTER – 4

SELF-ASSESSMENT AUDIOMETER

4.1 INTRODUCTION

The self-assessment audiometer project endeavors to revolutionize the field of hearing screening by introducing a novel, user-friendly solution. The traditional methods of audiology often involve cumbersome procedures conducted by audiologists, which can be both time-consuming and costly for individuals seeking assessment. Recognizing these limitations, the project aims to develop an innovative audiometer that empowers individuals to monitor their hearing health independently, without the need for specialized professionals.





Fig 4.1 Self-assessment audiometer

4.2. PORTABLE AUDIOMETER UTILIZING RASPBERRY PI

Central to the project is the development of a portable audiometer powered by Raspberry Pi technology. The Raspberry Pi serves as the core computing platform, offering a versatile and cost-effective solution for audiological testing. Leveraging its computational capabilities, the audiometer can generate pure tones across various frequency ranges and deliver them to the patient via headphones. The portability of the device ensures flexibility in conducting assessments in diverse settings, from clinical environments to remote or home-based settings.

4.3. PORTABLE AUDIOMETER (RASPBERRY PI) VS MOBILE HEARING TEST

FEATURES	PORTABLE AUDIOMETER (RASPBERRY PI)	MOBILE HEARING TEST
PICTURE	 <p>Fig 4.3.1 PORTABLE AUDIOMETER</p>	 <p>Fig 4.3.2. MOBILE HEARING TEST</p>
CALIBRATION LEVEL	Utilizes precise calibration of sound pressure levels (SPL)	Varies widely, dependent on the quality and calibration of mobile devices
TONE USED	Pure tones of varying frequencies	Pure tones or speech signals depending on the app
BENEFITS	Provides accurate and reliable results, customizable, integrates with various software	Highly accessible, convenient, and easy to use
CROSS-EAR LEVELS	Allows for detailed measurement of interaural attenuation (cross-hearing)	Limited capability to assess cross-ear hearing accurately
OPERATION	Requires basic technical knowledge for setup and operation	User-friendly apps, minimal technical knowledge required
PORTABILITY	Highly portable, but requires several components (Raspberry Pi, headphones, etc.)	Extremely portable, just requires a smartphone and headphones

COST	Moderate initial cost for hardware components	Low cost, often just the price of the app or free
EASE OF USE	Requires initial setup and familiarization	Intuitive interfaces, minimal setup required
FLEXIBILITY	High, can be customized for various audiometric tests	Moderate, depends on app features
DATA STORAGE	Saves results in CSV files, easy to manage and analyze	Saves results within the app or cloud storage, varies by app
TECHNOLOGY	Uses Raspberry Pi 3 B+, Python programming for automation	Utilizes smartphone technology and app-based interfaces

Table 2: Portable Audiometer (Raspberry Pi) Vs Mobile Hearing Test

4.4. Self-Assessment Audiometer For Testing: Possible And Not Possible Scenarios

1. Possible to Take Test:

- **Routine Home Screening:** Individuals can perform routine checks at home.
- **Convenient Monitoring:** Patients can monitor their hearing conveniently.
- **Remote Accessibility:** Accessible in remote areas without audiologists.
- **Cost-Effective Solution:** Affordable alternative to traditional tests.
- **Follow-Up Testing:** Easily conduct follow-up tests at home

2. Not Possible to Take Test:

- **Severe Hearing Impairments:** Profoundly deaf individuals may not respond accurately.
- **Pediatric Testing:** Young children and infants may not understand the test.

- **Complex Disorders:** Detailed diagnostic procedures are required for complex conditions.
- **Calibration and Accuracy:** Potential calibration issues affecting accuracy.
- **Technological Proficiency:** Difficulty for users lacking technological skills.
- **Environmental Noise:** Testing in noisy environments can lead to inaccurate results.
- **Professional Interpretation:** Some results may require professional interpretation.

4.5. EXISTING VS PROPOSED SELF-ASSESSMENT AUDIOMETERS

FEATURE	PORTABLE SELF-ASSESSMENT AUDIOMETER	TRADITIONAL AUDIOMETER
PICTURE	 Fig 4.5.1 PORTABLE AUDIOMETER	 Fig 4.5.2. AUDIOMETER
CALIBRATION LEVEL	Calibrated to a specific sound pressure level (SPL)	Requires professional calibration to maintain accuracy
TONE USED	Pure tones across various frequencies	Pure tones across various frequencies
BENEFITS	Cost-effective, convenient, and user-friendly	Highly accurate with professional oversight

CROSS-EAR LEVELS	Managed by software algorithms to minimize cross-hearing	Audiologist adjusts to prevent cross-hearing
OPERATION	Automated testing procedure with patient-controlled responses	Audiologist conducts and adjusts the test manually
PORTABILITY	Highly portable; consists of Raspberry Pi, headphones, button	Generally stationary, housed in clinics or hospitals
COST	Lower cost due to use of affordable components and open-source software	Higher cost due to sophisticated equipment and professional fees
ACCESSIBILITY	Accessible at home or in remote locations	Accessible primarily in medical or audiology facilities
EASE OF USE	Designed for ease of use by individuals with minimal training	Requires an audiologist to operate
FLEXIBILITY	Flexible; can be used in various settings and by different users	Limited to clinical settings and use by trained personnel
DATA STORAGE	Results saved automatically as CSV files with date/time stamps	Data typically recorded manually or through clinic software
TECHNOLOGY	Uses Raspberry Pi 3 B+, Python for software implementation	Utilizes advanced audiometric equipment and specialized software

TABLE 3: EXISTING VS PROPOSED SELF-ASSESSMENT AUDIOMETERS

CHAPTER – 5

PURE TONE AUDIOMETER

5.1. INTRODUCTION

A pure tone audiometer serves as a specialized tool for evaluating an individual's auditory capabilities. It administers pure tones of varying frequencies and strengths to each ear independently, enabling audiologists to ascertain the faintest sounds detectable by the person across different frequencies. This data is then utilized to generate an audiogram, a visual representation illustrating the individual's hearing thresholds across various frequencies.



Fig 5.1. Pure tone audiometer

5.2. PURPOSE OF PURE TONE AUDIOMETER

1. Assess hearing abilities accurately.
2. Diagnose various types and degrees of hearing loss.
3. Plan appropriate treatments, such as hearing aids or surgery.
4. Monitor changes in hearing sensitivity over time.
5. Screen for occupational and clinical hearing health.
6. Contribute to research and education in audiology.

5.3. TYPES OF TESTING



Fig 5.3. Types of testing

Pure tone audiometry comprises various testing modalities to evaluate different aspects of hearing. Here are the primary types:

1. Air Conduction Test:

This form of pure tone audiometry, widely used, entails delivering pure tones via headphones or speakers to determine the lowest audible sounds across a spectrum of frequencies. It examines the entire auditory pathway, from the external ear to the brain's processing of auditory stimuli.

2. Bone Conduction Test :

This testing method involves placing a bone vibrator behind the ear to directly stimulate the inner ear, bypassing the outer and middle ear. It helps discern whether hearing loss stems from issues in the outer or middle ear (conductive hearing loss) or the inner ear or auditory nerve (sensorineural hearing loss).

3. Speech Audiometry:

This evaluation gauges an individual's capacity to perceive and comprehend speech. It determines the minimum sound intensity at which an individual can recognize and repeat words, offering insights into speech perception and discrimination abilities.

4. Masking:

Masking is employed when one ear might inadvertently perceive sounds intended for the other ear, particularly in cases of significant hearing asymmetry. It entails presenting noise to the non-test ear to ensure that only the test ear responds to the pure tones.

5. Specialized Testing for Pediatric Population

Pediatric audiology may incorporate interactive behavioral testing methods, such as visual reinforcement audiometry (VRA) or conditioned play audiometry (CPA), which are more engaging for young children and yield reliable hearing thresholds.

6. Tympanometry:

Although not a pure tone test per se, tympanometry is often conducted alongside pure tone audiometry. It evaluates the mobility of the eardrum and the functionality of the middle ear, aiding in the diagnosis of conditions like otitis media or eustachian tube dysfunction.

5.4. AUDIOGRAM




Fig 5.4. Audiogram

An audiogram serves as a visual representation depicting an individual's hearing capacity across different frequencies and levels of intensity. It holds significant importance in audiology for evaluating and diagnosing hearing impairments. Typically, the graph illustrates frequencies ranging from low to high along the horizontal axis (x-axis), measured in Hertz (Hz), and sound intensity or volume along the vertical axis (y-axis), measured in decibels (dB).

The results of audiograms are derived from a series of tests wherein the individual responds to various tones or sounds. The audiologist documents the faintest sound audible to the person at each frequency, forming a profile of hearing thresholds. Offering insights into the type, severity, and pattern of hearing loss, the audiogram assists audiologists in devising appropriate interventions, such as recommending hearing aids or other necessary measures.



AUDIOPGRAM - All tests carried out in accordance with BSA protocols (Ref:BSA 2012)			
DATE OF TESTS: ____ / ____ / 2015	TESTED BY:	AUDIOLIST SIGNATURE:	
PATIENT NAME:	DATE OF BIRTH: ____ / ____ / ____	AGE:	PATIENT SIGNATURE:
AUDIOMETER:	SERIAL NO.:	ROOM SOUND LEVELS ACCEPTABLE? [] Yes [] No ____ db	
HEADPHONES:	CALIBRATED IN:	 YOUR WORLD NIHL AUDIOLOGY Accurate Diagnostic Audiometry	
PURE TONE AUDIOGRAM - RIGHT		PURE TONE AUDIOGRAM - LEFT	
			
AGE ADJUSTED SCORE: _____			
OTOSCOPY ACCEPTABLE [if No, see reverse] <input type="checkbox"/> Yes <input type="checkbox"/> No	EXPOSURE TO NOISE IN LAST 24 HOURS <input type="checkbox"/> Yes <input type="checkbox"/> No	SURGERY <input type="checkbox"/> Yes <input type="checkbox"/> No	
HEARING AID USER <input type="checkbox"/> Yes <input type="checkbox"/> No	PTA - REPEATABLE AND RELIABLE? <input type="checkbox"/> Yes <input type="checkbox"/> No	TINNITUS (Pulse/Warble) <input type="checkbox"/> Yes <input type="checkbox"/> No	

Fig 5.4.1. Audiogram Sheet

5.5. AUDIOGRAM CONFIGURATION

The structure of an audiogram chart illustrates an individual's hearing thresholds across various frequencies. Typical patterns encompass uniform (flat) loss, gradual (sloping) decline with higher frequencies predominantly affected, increasing (rising) loss with lower frequencies predominantly impacted, mid-frequency emphasis (U-shaped), and severe loss across all frequencies (corner).




Fig 5.5. Audiogram Configuration

5.6. LIMITATIONS

1. **Speech Understanding:** Audiograms don't measure how well someone understands speech, which can be affected by factors beyond hearing sensitivity.
2. **Sound Quality:** They don't capture sound quality or distortion, which can impact speech comprehension even if sounds are detected.
3. **Frequency Range:** Audiograms may not detect very low or very high-frequency losses, which can affect perception in specific situations.
4. **Localization:** They don't assess the ability to locate sounds, important for spatial awareness and safety.

5.7. APPLICATION

Pure tone audiometers are essential tools for:

1. **Hearing Assessment:** They accurately assess hearing sensitivity, aiding in diagnosis and classification of hearing loss.
2. **Treatment Monitoring:** Audiometers help monitor treatment effectiveness, such as hearing aid usage.
3. **Screening:** They're used for occupational and clinical hearing screenings.
4. **Research and Education:** Audiometers support research and education in audiology.

5.8. FEATURES

1. **Wide Frequency Range:** Covering low to high frequencies.
2. **Intensity Control:** Allows precise sound level adjustment.
3. **Ear Configurations:** Accommodates various testing methods, including air and bone conduction.
4. **Masking Capabilities:** Ensures accurate results, especially for asymmetric hearing loss.
5. **Data Management:** Offers storage and connectivity options for patient records.
6. **Portability:** Some models are portable for testing flexibility.

CHAPTER - 6

IMPLEMENTATION OF HARDWARE MODULES

6.1 INTRODUCTION

The portable pure tone audiometer utilizes a Raspberry Pi 3 B+, alongside a mouse and headphones. The Raspberry Pi operates on the Raspbian OS and is equipped with a 64-bit 1.2GHz quad-core processor and essential connectivity ports. Programmed in Python, it facilitates straightforward script development. Leveraging the Raspberry Pi's ALSA sound card, the audiometer produces pure tones spanning frequencies from 125 Hz to 8 kHz.

These tones are conveyed through a 3.5mm audio jack to HP102BK headphones worn by the patient. Testing initiates automatically with the right ear and is masked to prevent erroneous responses. The patient indicates tone perception by pressing the mouse button, prompting frequency and intensity adjustments through Python scripting.

6.2. BLOCK DIAGRAM

This initiative empowers individuals experiencing hearing challenges to self-assess their auditory abilities and ascertain the extent of their hearing range. The schematic representation of the envisaged portable audiometer is illustrated in Figure 6.2.




Fig 6.2. Block Diagram

COMPONENT DESCRIPTION

6.3.1 RASPBERRY PI

The Raspberry Pi, a line of compact single-board computers (SBCs), was developed by the Raspberry Pi Foundation in partnership with Broadcom in the United Kingdom. Initially designed to enrich fundamental computer science education in educational institutions, Raspberry Pi quickly expanded beyond its educational role, attracting a broad spectrum of users from enthusiasts to industrial pioneers.

Its cost-effectiveness, along with its modular structure and inclusion of standard interfaces like HDMI and USB, elevated Raspberry Pi into a flexible platform suitable for various uses, such as multimedia centers, robotics, and home automation setups.

The establishment of Raspberry Pi (Trading) Ltd, led by visionary Eben Upton, marked a pivotal moment in the evolution of the Raspberry Pi ecosystem. While maintaining its commitment to education, the Foundation pivoted towards commercial endeavors, ensuring the widespread availability of Raspberry Pi devices through manufacturing facilities in Wales, China, and Japan.

Surpassing the iconic ZX Spectrum in 2015, Raspberry Pi emerged as the best-selling British computer, solidifying its status as a beacon of innovation and accessibility in the computing realm. In 2021, the rebranding of Raspberry Pi (Trading) Ltd to Raspberry Pi Ltd reaffirmed its enduring mission to foster innovation and empower learners worldwide.



Fig 6.3.1 Raspberry Pi

6.3.1.1. FEATURES OF RASPBERRY PI:

The general features of this board are as follows:

1. **Affordability:** Raspberry Pi offers a low-cost computing solution.
2. **Versatility:** It suits a wide range of projects, from simple electronics to complex IoT applications.

3. **GPIO Pins:** GPIO pins enable interaction with external components.
4. **Community Support:** A large and active community provides ample resources and assistance.
5. **Education:** Raspberry Pi promotes learning programming and electronics, originally designed for educational purposes.




Fig 6.3.2 Pin Description of Raspberry Pi

6.3.1.2. RASPBERRY PI GPIO'S PIN

- ❖ GPIO (General Purpose Input/Output) pins enable interaction with external components. They act as programmable switches for controlling (output) or reading from (input) devices.
- ❖ These pins are grouped into sets for specific purposes like power, ground, digital I/O, and analog input.
- ❖ Raspberry Pi models offer varying numbers and configurations of GPIO pins, typically in a row of 40 (some models feature 26).
- ❖ GPIO pins facilitate a wide range of projects by allowing users to interface with sensors, LEDs, motors, and other peripherals.




Fig 6.3.3 Raspberry Pi GPIO's Pin Layout

6.3.1.3. RASPBERRY PI PIN DETAILS

1. **Power Pins:** Supply power to external devices (5V and 3.3V).
2. **Ground Pins:** Provide reference voltage for circuits.
3. **Digital I/O Pins:** Read inputs or control outputs.
4. **Analog Input Pins:** Read analog signals from sensors.
5. **Special Function Pins:** Handle specific communication protocols like I2C, SPI, UART, or PWM.




Fig 6.3.4 Raspberry Pi GPIO

6.4. RASPBERRY PI CASE:

A protective enclosure designed to house and shield the Raspberry Pi 3 B+ from physical damage and environmental hazards.

- **Features:**

1. Secure enclosure to safeguard the Raspberry Pi components.
2. Accessible openings for ports, GPIO pins, and camera/display connectors.

- **Specifications:**

- Material: Plastic or aluminum
- Design: Ventilation holes for heat dissipation, mounting holes for stability.



Fig 6.4. Raspberry Pi Case

6.5. MOUSE

The mouse is an input device used for interacting with graphical user interfaces on computers.

- **Features:**

1. Ergonomic design for comfortable use.
2. Responsive optical or laser sensor for precise cursor movement.

- **Specifications:**

- Connectivity: Wired or wireless (USB or Bluetooth)
- Buttons: Left, right, scroll wheel



Fig 6.5 Mouse

6.6. HEADPHONE:

The headphone delivers audio output to the user for listening to sound from the Raspberry Pi.

1. Features:

1. High-quality audio reproduction for immersive sound experience.
2. Comfortable design for extended wear.

2. Specifications:

- Type: Over-ear or in-ear
- Connectivity: Wired (3.5mm audio jack) or wireless (Bluetooth)
- Frequency Response: Typically ranges from 20Hz to 20kHz



Fig 6.6. Headphoney

6.6.1. TDH 49 HEADPHONE CHARACTERISTICS AND USAGE:

1. Tone Used:

- ❖ Pure tones across various frequencies for audiometric testing.

2. Positive Characteristics:

- ❖ Reliable and consistent in delivering precise sound stimuli.
- ❖ Closed design attenuates ambient noise for accurate testing.

3. Negative Characteristics:

- ❖ Bulky size and weight may cause discomfort during prolonged use.
- ❖ Closed design can lead to heat and moisture buildup, causing discomfort.

4. Cross Ear Level:

- ❖ Designed to minimize cross-hearing with good acoustic isolation.
- ❖ Proper positioning crucial for minimal sound leakage and accurate results.

5. Crossover:

- ❖ Potential feature for consistent sound quality across frequencies.
- ❖ Ensures accurate and reliable audiometric testing results.

6.7. ETHERNET CABLE:

The Ethernet cable establishes a wired network connection between devices for high-speed data transfer.

1. Features:

- Secure and stable network connectivity.
- Twisted pair construction for reduced interference.

2. Specifications:

- Cable Category: Cat5e, Cat6, Cat6a
- Length: Various lengths available
- Speed: Typically up to 1Gbps



Fig 6.7. Ethernet Cable:



Fig 6.8 Power Adapter

6.8. Power Supply

The power source delivers essential electrical energy to the Raspberry Pi and accompanying components of the audiometer, ensuring steady and dependable functionality of the system.

1. Characteristics:

- Provides consistent output voltage for reliable performance.
- Equipped with safety measures like overload and short-circuit protection.

2. Specifications:

- Input: AC 100-240V, 50/60Hz
- Output: DC 5V, 2.5A
- Connector Type: Micro USB

CHAPTER – 7

SOFTWARE IMPLEMENTATION

7.1 INTRODUCTION

In Technology, the word implementation of software is used to describe the software tools implemented in this project.

7.2. RASPBERRY PI OS (RASPBIAN)

Installing Raspberry Pi OS involves preparing a microSD card and setting up your Raspberry Pi. Here's a quick overview

- ❖ Installing Raspberry Pi OS via Raspberry Pi Imager offers a swift and straightforward method to set up the operating system on a microSD card, making it compatible with your Raspberry Pi for immediate use.

https://downloads.raspberrypi.org/imager/imager_latest.exe




Fig 7.2 Installing Raspbian.

1. **Download Raspberry Pi Imager:** This tool writes the OS to your SD card. Get it from the Raspberry Pi website.
2. **Prepare microSD card:** Use a minimum 8GB class 10 card. The imager will erase it during setup.
3. **Flash OS to SD card:** Use the imager to select Raspberry Pi OS (formerly Raspbian) and your SD card. The imager will write the OS to the card.

4. **Setup Raspberry Pi:** Insert the SD card, connect a monitor, keyboard, mouse, and power supply to your Pi. You can also connect to Wi-Fi during initial setup.




Fig 7.2.1 The Raspberry Pi desktop.

6.3. VNC: Remote access a Raspberry Pi

To install VNC Viewer on your computer and connect to your Raspberry Pi:

1. **Download VNC Viewer:** Visit the RealVNC website and download the VNC Viewer application for your operating system.
2. **Install VNC Viewer:** Follow the on-screen instructions to install VNC Viewer on your computer.
3. **Configure Raspberry Pi:** On your Raspberry Pi, enable VNC Server. You can do this by opening a terminal and running `sudo raspi-config`, navigating to **Interfacing Options**, and enabling VNC.
4. **Find Raspberry Pi IP Address:** Use `ifconfig` on the Raspberry Pi terminal to find its IP address.
5. **Connect with VNC Viewer:** Open VNC Viewer on your computer and enter the Raspberry Pi's IP address. Click connect and enter the username and password for your Raspberry Pi when prompted.
6. **Access Raspberry Pi Desktop:** You should now be connected to your Raspberry Pi's desktop environment using VNC Viewer.




Fig 7.3 VNC: Remote Access to a Raspberry Pi

7.4. Python on Raspberry Pi: Thonny Editor

The Raspberry Pi OS includes a basic Python editor called "Thonny", providing a straightforward environment for writing and executing Python code directly on the Raspberry Pi.




Fig 7.4 Thonny ide

CHAPTER – 8

VALIDATION AND FEEDBACK: MEETING WITH ENT DOCTOR

Consultation with Dr. B. Harsha Vardhan, M.S. ENT (Otolaryngology)

8.1. INTRODUCTION

To develop a portable audiometer for self-assessment of hearing loss, we sought input from Dr. B. Harsha Vardhan, an ENT expert, to gather practical insights and evaluate the device's effectiveness. Dr. Vardhan provided valuable guidance and feedback on the project.

8.2. CONSULTATION WITH DR. B. HARSHA VARDHAN, M.S. ENT (OTOLARYNGOLOGY)




Fig 8.2. Consulting with the Doctor

- ◆ We met Dr. B. Harsha Vardhan, an ENT specialist based in Mahalingapuram, Pollachi.
- ◆ During the meeting, we discussed the project and demonstrated the functionality of the pure tone audiometer.
- ◆ Dr. Vardhan praised the project, considering it highly promising.
- ◆ He expressed interest in purchasing the pure tone audiometer once it is completed.

8.3. Recommendation to Meet an Audiologist:

- ◆ Dr. Vardhan recommended meeting with an audiologist in Pollachi for further validation of the device.
- ◆ He provided contact details of an audiologist in Pollachi for additional consultation.

8.4. Discussion with the ENT Doctor:

- ◆ We spoke with the audiologist and discussed the project in detail.
- ◆ We learned that the market value of a commercial pure tone audiometer is approximately 1.50 lakhs.
- ◆ The audiologist advised us to test the self-assembled pure tone audiometer with individuals who are hearing impaired for further validation.
- ◆ Feedback was given on making minor adjustments to the audiogram's x and y-axis for improved accuracy.

8.5. Testing the Self-Assembled Pure Tone Audiometer:

- ◆ We met with an audiologist in Pollachi who had access to a commercial pure tone audiometer.
- ◆ Comparative tests were conducted between the self-assembled audiometer and the commercial device.
- ◆ The results indicated that the self-assembled device performed well, with only slight variations in audio tones

8.6. Feedback and Conclusion:

- ◆ Dr. Vardhan appreciated the innovative project and participated in testing the audiometer.
- ◆ The project received positive feedback overall, with suggestions for minor adjustments to enhance accuracy.
- ◆ Collaboration with medical professionals validated the effectiveness of the self-assembled pure tone audiometer, reinforcing its potential utility.

CHAPTER - 9

MEETING WITH AUDIOLOGIST BALACHANDRAN, ASLP

Consultation with Dr. P. Balachandran, ASLP

9.1. INTRODUCTION:

During a recent visit to Pollachi, I had the opportunity to meet with Audiologist Balachandran, ASLP. We engaged in an in-depth discussion about audiometric equipment and testing procedures for over an hour.

9.2. MEETING WITH AUDIOLOGIST BALACHANDRAN ASLP IN POLLACHI

In Pollachi, I met with Audiologist Balachandran, ASLP, for an enlightening conversation lasting over an hour. He introduced me to essential audiology testing equipment, including the Auditivio Audiometer (Endeavour), an Indian product priced at ₹50,000, and the imported Flute Inventis Tympanometer, which costs ₹4.50 lakhs.




Fig 9.2. Meeting with Audiologist

9.3. PROJECT DISCUSSION AND RECOMMENDATIONS

Balachandran provided valuable insights and suggestions for enhancing our project. He emphasized the importance of including a bone conduction test to accurately plot the audiogram and describe hearing loss. He also recommended replacing the standard headphones with TDH49, which are specifically designed for audiometry tests. Furthermore, he suggested printing the audiogram sheet for proper documentation.



Fig 9.3. Audiometry Test

9.4. QUESTIONS AND AUDIOMETER INSIGHTS

During our conversation, Balachandran posed several insightful questions regarding audiometry:

- ◆ He inquired about the differences between our project and a pure tone audiometer.
- ◆ He discussed the disparity between a hearing test app and a self-assessment audiometry test.
- ◆ He questioned the significance of intensity changes in our project.
- ◆ He suggested studying the characteristics of the TDH49 headphones, including their crossover and ear levels.

- ◆ He prompted us to explore calibration levels for sound pressure levels in our self-assessment portable audiometer.
- ◆ He compared the cost-effectiveness of our project with a pure tone audiometer.
- ◆ He raised queries about the intensity and time delay in our project.
- ◆ He advised researching the maximum frequency intensity of an audiometer.

9.5. GUIDANCE AND LEARNING

Balachandran recommended the book "Auditory Diagnosis" by Robert and requested a hard copy of our project documentation. We provided him with detailed study material. He also explained the Mackenzie audiometry test, highlighting its significance.

9.6. KEY FEEDBACK FROM AUDIOLOGIST BALACHANDRAN

- ◆ **Bone Conduction Tests:** Essential for accurate audiogram plotting and describing hearing loss.
- ◆ **TDH49 Headphones:** Recommended replacing standard headphones with TDH49 for audiometry tests.
- ◆ **Printed Audiogram Sheets:** Necessary for proper documentation.
- ◆ **Calibration Levels:** Explore and implement sound pressure level calibration for accuracy.
- ◆ **Intensity and Time Delay Adjustments:** Needed for better test accuracy.
- ◆ **Maximum Frequency Intensity:** Research the maximum frequency intensity levels.
- ◆ **TDH49 Headphone Characteristics:** Study crossover and ear level responses.
- ◆ **Cost-Effectiveness:** Compare the project's cost-effectiveness with pure tone audiometers.

- ◆ **Mackenzie Audiometry Test:** Consider incorporating or adapting principles from this test
 - ◆ **High-Frequency Audiometry Tests:** Extend capabilities to include high-frequency audiometry tests.
 - ◆ **Detailed Documentation and User Manual:** Develop comprehensive documentation and a user manual.
 - ◆ **Enhanced User Interface:** Improve user interface for better usability.
 - ◆ **Real-Time Data Analysis and Feedback:** Integrate real-time data analysis and feedback features for immediate results.



Fig 9.6. Auditivo Audiometer

9.7. EXPLORATION AND FUTURE SCOPE

Balachandran educated us on high-frequency audiometry tests and provided guidance for further development of our project. He offered detailed insights into masking tests, cross-level intensity, and various types of audiometry, and explained the concept of hearing loss.

9.8. CONCLUSION

The meeting with Balachandran was incredibly enlightening. We gained valuable insights into potential improvements and future prospects for our project. His expertise in audiology testing and his recommendations were instrumental in enhancing our understanding and guiding our project's direction.

CHAPTER – 10

RESULTS

10.1 INTRODUCTION

This chapter explains the testing of the working of hardware, software, raspberry pi (Pure Tone Audiometry) of the project.

10.2. TESTING OF FINAL WORKING MODEL:

The integration of the audiogram for assessing hearing loss was effectively achieved within the Python framework, aided by a computer, a Raspberry Pi, and headphones. Utilizing the sound card of the Raspberry Pi 3 B+, pure tones of different frequencies and strengths were produced. These tones were subsequently delivered to the headphones worn by the hearing-impaired individual. The individual's reactions, triggered by clicking the left mouse button, were captured and logged by the Python program through the Raspberry Pi configuration.

Before commencing the audiometric test, individuals with hearing impairment were provided with essential software instructions:

1. Functional Verification:

The portable audiometer effectively conducted pure tone audiometry tests, generating automatic audiograms based on patient responses.

2. Data Storage and Analysis:

Test results were stored in CSV files, enabling detailed analysis of patient responses and comparison with traditional audiometry results.

3. Visualizations:

Audiogram charts provided visual representations of hearing levels across different frequencies, aiding interpretation of test results for both patients and healthcare professionals.

4. Accuracy and Reliability:

Validation testing demonstrated the audiometer's accuracy and reliability, with results comparable to those obtained from conventional audiometers.

5. Feedback Incorporation:

Feedback from medical professionals informed refinements to the audiometer's performance, enhancing its usability and effectiveness in clinical settings.

10.3 WORKING MODEL

1. Components:

Raspberry Pi 3 B+, computer, patient response button, and headphones.

2. Sound Signal Delivery:

Headphones deliver sound signals of varying frequencies and volumes.

3. Patient Response:

Patients respond using a mouse button, providing feedback for recording.

4. Data Processing:

Real-time processing captures frequency, volume, and reaction time data.

5. Output:

Generates audiograms and visualizations for easy interpretation of hearing levels.




Fig 10.3 Overall Initialization Pure Tone Audiometer

10.4. RESULT ANALYSIS:

Using the data stored in the CSV file, the audiogram was produced to illustrate the hearing thresholds for both the right and left ears. The graph displayed the frequency, measured in Hertz, along the X-axis, and the sound intensity, measured in decibels (dB), along the Y-axis.

Fig 10.4.1. Python Code in Thonny Ide

```

Thonny - C:\Users\moham\Downloads\Mini_Project(Batch_08)\Mini Project Sample test\43(perfect_instruction).py @ 239 : 1
File Edit View Run Tools Help
Instructions
43(perfect
1
2
3
4
5
6
7
8
9
10
11
12
13 # Use interactive backend for displaying plots in a separate window

Portable Self Assessment Audiometer

1. Put on your headphones.
2. Click the left mouse button when you hear pulsing sounds.
3. The test will run for the right ear, playing sounds at different frequencies and volumes.
4. After each sound, click when you hear it.
5. Once finished, visualizations and files will be available for review.

Start Test

```

Shell >

```
>>> %Run '43(perfect_instruction).py'
```

Fig 10.4.2. Instruction of the project

```

Thonny - C:\Users\moham\Downloads\Mini_Project(Batch_08)\Mini Project Sample test\43(perfect_instruction).py @ 13 : 47
File Edit View Run Tools Help
43(perfect_instruction).py
1 import argparse
2 from datetime import datetime, timedelta
Shell >
>>> %Run '43(perfect_instruction).py'
Testing right ear...
Playing frequency: 125 Hz at volume: 0 dB for right ear
Recording event: [125, 0, datetime.datetime(2024, 5, 9, 18, 18, 52, 767101), datetime.datetime(2024, 5, 9, 18, 18, 52, 988576)]
Recording event: [125, 0, datetime.datetime(2024, 5, 9, 18, 18, 52, 767101), datetime.datetime(2024, 5, 9, 18, 18, 54, 377172)]
Playing frequency: 250 Hz at volume: 0 dB for right ear
Recording event: [250, 0, datetime.datetime(2024, 5, 9, 18, 18, 57, 218996), datetime.datetime(2024, 5, 9, 18, 18, 58, 565294)]
Playing frequency: 500 Hz at volume: 0 dB for right ear
Playing frequency: 500 Hz at volume: 10 dB for right ear
Recording event: [500, 10, datetime.datetime(2024, 5, 9, 18, 19, 4, 93461), datetime.datetime(2024, 5, 9, 18, 19, 4, 737185)]
Playing frequency: 1000 Hz at volume: 0 dB for right ear
Playing frequency: 1000 Hz at volume: 10 dB for right ear
Playing frequency: 1000 Hz at volume: 20 dB for right ear
Recording event: [1000, 20, datetime.datetime(2024, 5, 9, 18, 19, 13, 411136), datetime.datetime(2024, 5, 9, 18, 19, 14, 225200)]
Playing frequency: 2000 Hz at volume: 0 dB for right ear
Playing frequency: 2000 Hz at volume: 10 dB for right ear
Recording event: [2000, 10, datetime.datetime(2024, 5, 9, 18, 19, 20, 281548), datetime.datetime(2024, 5, 9, 18, 19, 21, 619935)]

```

Fig 10.4.3. Patient Response of the Audiometer

Thonny - C:\Users\moham\Downloads\Mini_Project\Batch_08\Mini Project Sample test\43(perfect_instruction).py @ 13 : 1

File Edit View Run Tools Help

Program arguments:

```
43(perfect_instruction).py
```

```

1 import argparse
2 from datetime import datetime
3 from time import sleep
4 import numbers
5 import pyaudio
6 import threading
7 import pandas as pd
8 import matplotlib.pyplot as plt
9 from pydub import AudioSegment
10 from pydub.playback import play
11 import tkinter as tk
12
13 # Use integer base 10 for volume
14
15 def main():
16     # Set up audio recording
17     p = pyaudio.PyAudio()
18     stream = p.open(format=pyaudio.paInt16, channels=1, rate=44100, input=True, frames_per_buffer=1024)
19
20     # Set up plotting
21     fig, ax = plt.subplots()
22     ax.set_title("Audiogram for right ear")
23     ax.set_xlabel("Pitch (frequency in Hz)")
24     ax.set_ylabel("Hearing Level in decibels (volume in dB)")
25     ax.set_xlim(0, 8000)
26     ax.set_ylim(-10, 90)
27     ax.grid(True)
28
29     # Set up color mapping for hearing loss levels
30     colors = {
31         "Normal Hearing (0-15 dB)": "#2ca02c",
32         "Slight Hearing Loss (16-25 dB)": "#ffff99",
33         "Mild Hearing Loss (26-40 dB)": "#ffcc70",
34         "Moderate Hearing Loss (41-55 dB)": "#ff9999",
35         "Moderately Severe Hearing Loss (56-70 dB)": "#9970c0",
36         "Severe Hearing Loss (71-90 dB)": "#c07070",
37         "Profound Hearing Loss (91 dB or greater)": "#707070"
38     }
39
40     # Plot the audiogram
41     while True:
42         data = stream.read(1024)
43         play(data)
44         spectrum = np.abs(np.fft.rfft(data))
45         threshold = np.percentile(spectrum, 90)
46         threshold_db = 20 * np.log10(threshold)
47         ax.plot([0, 8000], [threshold_db, threshold_db], 'k')
48         ax.plot([0, 8000], [threshold_db + 10, threshold_db + 10], 'k')
49         ax.plot([0, 8000], [threshold_db - 10, threshold_db - 10], 'k')
50
51         # Create a heatmap of the audiogram
52         heatmap = np.zeros((8000, 100), dtype=np.uint8)
53         for f in range(0, 8000):
54             for t in range(0, 100):
55                 if f < 1000:
56                     if t > threshold_db:
57                         heatmap[f, t] = 1
58                 else:
59                     heatmap[f, t] = 0
60             else:
61                 if f < 2000:
62                     if t > threshold_db + 10:
63                         heatmap[f, t] = 1
64                     else:
65                         heatmap[f, t] = 0
66                 else:
67                     if t > threshold_db - 10:
68                         heatmap[f, t] = 1
69                     else:
70                         heatmap[f, t] = 0
71
72         # Create a heatmap of the audiogram
73         heatmap = np.zeros((8000, 100), dtype=np.uint8)
74         for f in range(0, 8000):
75             for t in range(0, 100):
76                 if f < 1000:
77                     if t > threshold_db:
78                         heatmap[f, t] = 1
79                 else:
80                     heatmap[f, t] = 0
81             else:
82                 if f < 2000:
83                     if t > threshold_db + 10:
84                         heatmap[f, t] = 1
85                     else:
86                         heatmap[f, t] = 0
87                 else:
88                     if t > threshold_db - 10:
89                         heatmap[f, t] = 1
90                     else:
91                         heatmap[f, t] = 0
92
93         # Create a heatmap of the audiogram
94         heatmap = np.zeros((8000, 100), dtype=np.uint8)
95         for f in range(0, 8000):
96             for t in range(0, 100):
97                 if f < 1000:
98                     if t > threshold_db:
99                         heatmap[f, t] = 1
100                else:
101                    heatmap[f, t] = 0
102            else:
103                if f < 2000:
104                    if t > threshold_db + 10:
105                        heatmap[f, t] = 1
106                    else:
107                        heatmap[f, t] = 0
108                else:
109                    if t > threshold_db - 10:
110                        heatmap[f, t] = 1
111                    else:
112                        heatmap[f, t] = 0
113
114         # Create a heatmap of the audiogram
115         heatmap = np.zeros((8000, 100), dtype=np.uint8)
116         for f in range(0, 8000):
117             for t in range(0, 100):
118                 if f < 1000:
119                     if t > threshold_db:
120                         heatmap[f, t] = 1
121                 else:
122                     heatmap[f, t] = 0
123             else:
124                 if f < 2000:
125                     if t > threshold_db + 10:
126                         heatmap[f, t] = 1
127                     else:
128                         heatmap[f, t] = 0
129                 else:
130                     if t > threshold_db - 10:
131                         heatmap[f, t] = 1
132                     else:
133                         heatmap[f, t] = 0
134
135         # Create a heatmap of the audiogram
136         heatmap = np.zeros((8000, 100), dtype=np.uint8)
137         for f in range(0, 8000):
138             for t in range(0, 100):
139                 if f < 1000:
140                     if t > threshold_db:
141                         heatmap[f, t] = 1
142                 else:
143                     heatmap[f, t] = 0
144             else:
145                 if f < 2000:
146                     if t > threshold_db + 10:
147                         heatmap[f, t] = 1
148                     else:
149                         heatmap[f, t] = 0
150                 else:
151                     if t > threshold_db - 10:
152                         heatmap[f, t] = 1
153                     else:
154                         heatmap[f, t] = 0
155
156         # Create a heatmap of the audiogram
157         heatmap = np.zeros((8000, 100), dtype=np.uint8)
158         for f in range(0, 8000):
159             for t in range(0, 100):
160                 if f < 1000:
161                     if t > threshold_db:
162                         heatmap[f, t] = 1
163                 else:
164                     heatmap[f, t] = 0
165             else:
166                 if f < 2000:
167                     if t > threshold_db + 10:
168                         heatmap[f, t] = 1
169                     else:
170                         heatmap[f, t] = 0
171                 else:
172                     if t > threshold_db - 10:
173                         heatmap[f, t] = 1
174                     else:
175                         heatmap[f, t] = 0
176
177         # Create a heatmap of the audiogram
178         heatmap = np.zeros((8000, 100), dtype=np.uint8)
179         for f in range(0, 8000):
180             for t in range(0, 100):
181                 if f < 1000:
182                     if t > threshold_db:
183                         heatmap[f, t] = 1
184                 else:
185                     heatmap[f, t] = 0
186             else:
187                 if f < 2000:
188                     if t > threshold_db + 10:
189                         heatmap[f, t] = 1
190                     else:
191                         heatmap[f, t] = 0
192                 else:
193                     if t > threshold_db - 10:
194                         heatmap[f, t] = 1
195                     else:
196                         heatmap[f, t] = 0
197
198         # Create a heatmap of the audiogram
199         heatmap = np.zeros((8000, 100), dtype=np.uint8)
200         for f in range(0, 8000):
201             for t in range(0, 100):
202                 if f < 1000:
203                     if t > threshold_db:
204                         heatmap[f, t] = 1
205                 else:
206                     heatmap[f, t] = 0
207             else:
208                 if f < 2000:
209                     if t > threshold_db + 10:
210                         heatmap[f, t] = 1
211                     else:
212                         heatmap[f, t] = 0
213                 else:
214                     if t > threshold_db - 10:
215                         heatmap[f, t] = 1
216                     else:
217                         heatmap[f, t] = 0
218
219         # Create a heatmap of the audiogram
220         heatmap = np.zeros((8000, 100), dtype=np.uint8)
221         for f in range(0, 8000):
222             for t in range(0, 100):
223                 if f < 1000:
224                     if t > threshold_db:
225                         heatmap[f, t] = 1
226                 else:
227                     heatmap[f, t] = 0
228             else:
229                 if f < 2000:
230                     if t > threshold_db + 10:
231                         heatmap[f, t] = 1
232                     else:
233                         heatmap[f, t] = 0
234                 else:
235                     if t > threshold_db - 10:
236                         heatmap[f, t] = 1
237                     else:
238                         heatmap[f, t] = 0
239
240         # Create a heatmap of the audiogram
241         heatmap = np.zeros((8000, 100), dtype=np.uint8)
242         for f in range(0, 8000):
243             for t in range(0, 100):
244                 if f < 1000:
245                     if t > threshold_db:
246                         heatmap[f, t] = 1
247                 else:
248                     heatmap[f, t] = 0
249             else:
250                 if f < 2000:
251                     if t > threshold_db + 10:
252                         heatmap[f, t] = 1
253                     else:
254                         heatmap[f, t] = 0
255                 else:
256                     if t > threshold_db - 10:
257                         heatmap[f, t] = 1
258                     else:
259                         heatmap[f, t] = 0
260
261         # Create a heatmap of the audiogram
262         heatmap = np.zeros((8000, 100), dtype=np.uint8)
263         for f in range(0, 8000):
264             for t in range(0, 100):
265                 if f < 1000:
266                     if t > threshold_db:
267                         heatmap[f, t] = 1
268                 else:
269                     heatmap[f, t] = 0
270             else:
271                 if f < 2000:
272                     if t > threshold_db + 10:
273                         heatmap[f, t] = 1
274                     else:
275                         heatmap[f, t] = 0
276                 else:
277                     if t > threshold_db - 10:
278                         heatmap[f, t] = 1
279                     else:
280                         heatmap[f, t] = 0
281
282         # Create a heatmap of the audiogram
283         heatmap = np.zeros((8000, 100), dtype=np.uint8)
284         for f in range(0, 8000):
285             for t in range(0, 100):
286                 if f < 1000:
287                     if t > threshold_db:
288                         heatmap[f, t] = 1
289                 else:
290                     heatmap[f, t] = 0
291             else:
292                 if f < 2000:
293                     if t > threshold_db + 10:
294                         heatmap[f, t] = 1
295                     else:
296                         heatmap[f, t] = 0
297                 else:
298                     if t > threshold_db - 10:
299                         heatmap[f, t] = 1
300                     else:
301                         heatmap[f, t] = 0
302
303         # Create a heatmap of the audiogram
304         heatmap = np.zeros((8000, 100), dtype=np.uint8)
305         for f in range(0, 8000):
306             for t in range(0, 100):
307                 if f < 1000:
308                     if t > threshold_db:
309                         heatmap[f, t] = 1
310                 else:
311                     heatmap[f, t] = 0
312             else:
313                 if f < 2000:
314                     if t > threshold_db + 10:
315                         heatmap[f, t] = 1
316                     else:
317                         heatmap[f, t] = 0
318                 else:
319                     if t > threshold_db - 10:
320                         heatmap[f, t] = 1
321                     else:
322                         heatmap[f, t] = 0
323
324         # Create a heatmap of the audiogram
325         heatmap = np.zeros((8000, 100), dtype=np.uint8)
326         for f in range(0, 8000):
327             for t in range(0, 100):
328                 if f < 1000:
329                     if t > threshold_db:
330                         heatmap[f, t] = 1
331                 else:
332                     heatmap[f, t] = 0
333             else:
334                 if f < 2000:
335                     if t > threshold_db + 10:
336                         heatmap[f, t] = 1
337                     else:
338                         heatmap[f, t] = 0
339                 else:
340                     if t > threshold_db - 10:
341                         heatmap[f, t] = 1
342                     else:
343                         heatmap[f, t] = 0
344
345         # Create a heatmap of the audiogram
346         heatmap = np.zeros((8000, 100), dtype=np.uint8)
347         for f in range(0, 8000):
348             for t in range(0, 100):
349                 if f < 1000:
350                     if t > threshold_db:
351                         heatmap[f, t] = 1
352                 else:
353                     heatmap[f, t] = 0
354             else:
355                 if f < 2000:
356                     if t > threshold_db + 10:
357                         heatmap[f, t] = 1
358                     else:
359                         heatmap[f, t] = 0
360                 else:
361                     if t > threshold_db - 10:
362                         heatmap[f, t] = 1
363                     else:
364                         heatmap[f, t] = 0
365
366         # Create a heatmap of the audiogram
367         heatmap = np.zeros((8000, 100), dtype=np.uint8)
368         for f in range(0, 8000):
369             for t in range(0, 100):
370                 if f < 1000:
371                     if t > threshold_db:
372                         heatmap[f, t] = 1
373                 else:
374                     heatmap[f, t] = 0
375             else:
376                 if f < 2000:
377                     if t > threshold_db + 10:
378                         heatmap[f, t] = 1
379                     else:
380                         heatmap[f, t] = 0
381                 else:
382                     if t > threshold_db - 10:
383                         heatmap[f, t] = 1
384                     else:
385                         heatmap[f, t] = 0
386
387         # Create a heatmap of the audiogram
388         heatmap = np.zeros((8000, 100), dtype=np.uint8)
389         for f in range(0, 8000):
390             for t in range(0, 100):
391                 if f < 1000:
392                     if t > threshold_db:
393                         heatmap[f, t] = 1
394                 else:
395                     heatmap[f, t] = 0
396             else:
397                 if f < 2000:
398                     if t > threshold_db + 10:
399                         heatmap[f, t] = 1
400                     else:
401                         heatmap[f, t] = 0
402                 else:
403                     if t > threshold_db - 10:
404                         heatmap[f, t] = 1
405                     else:
406                         heatmap[f, t] = 0
407
408         # Create a heatmap of the audiogram
409         heatmap = np.zeros((8000, 100), dtype=np.uint8)
410         for f in range(0, 8000):
411             for t in range(0, 100):
412                 if f < 1000:
413                     if t > threshold_db:
414                         heatmap[f, t] = 1
415                 else:
416                     heatmap[f, t] = 0
417             else:
418                 if f < 2000:
419                     if t > threshold_db + 10:
420                         heatmap[f, t] = 1
421                     else:
422                         heatmap[f, t] = 0
423                 else:
424                     if t > threshold_db - 10:
425                         heatmap[f, t] = 1
426                     else:
427                         heatmap[f, t] = 0
428
429         # Create a heatmap of the audiogram
430         heatmap = np.zeros((8000, 100), dtype=np.uint8)
431         for f in range(0, 8000):
432             for t in range(0, 100):
433                 if f < 1000:
434                     if t > threshold_db:
435                         heatmap[f, t] = 1
436                 else:
437                     heatmap[f, t] = 0
438             else:
439                 if f < 2000:
440                     if t > threshold_db + 10:
441                         heatmap[f, t] = 1
442                     else:
443                         heatmap[f, t] = 0
444                 else:
445                     if t > threshold_db - 10:
446                         heatmap[f, t] = 1
447                     else:
448                         heatmap[f, t] = 0
449
450         # Create a heatmap of the audiogram
451         heatmap = np.zeros((8000, 100), dtype=np.uint8)
452         for f in range(0, 8000):
453             for t in range(0, 100):
454                 if f < 1000:
455                     if t > threshold_db:
456                         heatmap[f, t] = 1
457                 else:
458                     heatmap[f, t] = 0
459             else:
460                 if f < 2000:
461                     if t > threshold_db + 10:
462                         heatmap[f, t] = 1
463                     else:
464                         heatmap[f, t] = 0
465                 else:
466                     if t > threshold_db - 10:
467                         heatmap[f, t] = 1
468                     else:
469                         heatmap[f, t] = 0
470
471         # Create a heatmap of the audiogram
472         heatmap = np.zeros((8000, 100), dtype=np.uint8)
473         for f in range(0, 8000):
474             for t in range(0, 100):
475                 if f < 1000:
476                     if t > threshold_db:
477                         heatmap[f, t] = 1
478                 else:
479                     heatmap[f, t] = 0
480             else:
481                 if f < 2000:
482                     if t > threshold_db + 10:
483                         heatmap[f, t] = 1
484                     else:
485                         heatmap[f, t] = 0
486                 else:
487                     if t > threshold_db - 10:
488                         heatmap[f, t] = 1
489                     else:
490                         heatmap[f, t] = 0
491
492         # Create a heatmap of the audiogram
493         heatmap = np.zeros((8000, 100), dtype=np.uint8)
494         for f in range(0, 8000):
495             for t in range(0, 100):
496                 if f < 1000:
497                     if t > threshold_db:
498                         heatmap[f, t] = 1
499                 else:
500                     heatmap[f, t] = 0
501             else:
502                 if f < 2000:
503                     if t > threshold_db + 10:
504                         heatmap[f, t] = 1
505                     else:
506                         heatmap[f, t] = 0
507                 else:
508                     if t > threshold_db - 10:
509                         heatmap[f, t] = 1
510                     else:
511                         heatmap[f, t] = 0
512
513         # Create a heatmap of the audiogram
514         heatmap = np.zeros((8000, 100), dtype=np.uint8)
515         for f in range(0, 8000):
516             for t in range(0, 100):
517                 if f < 1000:
518                     if t > threshold_db:
519                         heatmap[f, t] = 1
520                 else:
521                     heatmap[f, t] = 0
522             else:
523                 if f < 2000:
524                     if t > threshold_db + 10:
525                         heatmap[f, t] = 1
526                     else:
527                         heatmap[f, t] = 0
528                 else:
529                     if t > threshold_db - 10:
530                         heatmap[f, t] = 1
531                     else:
532                         heatmap[f, t] = 0
533
534         # Create a heatmap of the audiogram
535         heatmap = np.zeros((8000, 100), dtype=np.uint8)
536         for f in range(0, 8000):
537             for t in range(0, 100):
538                 if f < 1000:
539                     if t > threshold_db:
540                         heatmap[f, t] = 1
541                 else:
542                     heatmap[f, t] = 0
543             else:
544                 if f < 2000:
545                     if t > threshold_db + 10:
546                         heatmap[f, t] = 1
547                     else:
548                         heatmap[f, t] = 0
549                 else:
550                     if t > threshold_db - 10:
551                         heatmap[f, t] = 1
552                     else:
553                         heatmap[f, t] = 0
554
555         # Create a heatmap of the audiogram
556         heatmap = np.zeros((8000, 100), dtype=np.uint8)
557         for f in range(0, 8000):
558             for t in range(0, 100):
559                 if f < 1000:
560                     if t > threshold_db:
561                         heatmap[f, t] = 1
562                 else:
563                     heatmap[f, t] = 0
564             else:
565                 if f < 2000:
566                     if t > threshold_db + 10:
567                         heatmap[f, t] = 1
568                     else:
569                         heatmap[f, t] = 0
570                 else:
571                     if t > threshold_db - 10:
572                         heatmap[f, t] = 1
573                     else:
574                         heatmap[f, t] = 0
575
576         # Create a heatmap of the audiogram
577         heatmap = np.zeros((8000, 100), dtype=np.uint8)
578         for f in range(0, 8000):
579             for t in range(0, 100):
580                 if f < 1000:
581                     if t > threshold_db:
582                         heatmap[f, t] = 1
583                 else:
584                     heatmap[f, t] = 0
585             else:
586                 if f < 2000:
587                     if t > threshold_db + 10:
588                         heatmap[f, t] = 1
589                     else:
590                         heatmap[f, t] = 0
591                 else:
592                     if t > threshold_db - 10:
593                         heatmap[f, t] = 1
594                     else:
595                         heatmap[f, t] = 0
596
597         # Create a heatmap of the audiogram
598         heatmap = np.zeros((8000, 100), dtype=np.uint8)
599         for f in range(0, 8000):
600             for t in range(0, 100):
601                 if f < 1000:
602                     if t > threshold_db:
603                         heatmap[f, t] = 1
604                 else:
605                     heatmap[f, t] = 0
606             else:
607                 if f < 2000:
608                     if t > threshold_db + 10:
609                         heatmap[f, t] = 1
610                     else:
611                         heatmap[f, t] = 0
612                 else:
613                     if t > threshold_db - 10:
614                         heatmap[f, t] = 1
615                     else:
616                         heatmap[f, t] = 0
617
618         # Create a heatmap of the audiogram
619         heatmap = np.zeros((8000, 100), dtype=np.uint8)
620         for f in range(0, 8000):
621             for t in range(0, 100):
622                 if f < 1000:
623                     if t > threshold_db:
624                         heatmap[f, t] = 1
625                 else:
626                     heatmap[f, t] = 0
627             else:
628                 if f < 2000:
629                     if t > threshold_db + 10:
630                         heatmap[f, t] = 1
631                     else:
632                         heatmap[f, t] = 0
633                 else:
634                     if t > threshold_db - 10:
635                         heatmap[f, t] = 1
636                     else:
637                         heatmap[f, t] = 0
638
639         # Create a heatmap of the audiogram
640         heatmap = np.zeros((8000, 100), dtype=np.uint8)
641         for f in range(0, 8000):
642             for t in range(0, 100):
643                 if f < 1000:
644                     if t > threshold_db:
645                         heatmap[f, t] = 1
646                 else:
647                     heatmap[f, t] = 0
648             else:
649                 if f < 2000:
650                     if t > threshold_db + 10:
651                         heatmap[f, t] = 1
652                     else:
653                         heatmap[f, t] = 0
654                 else:
655                     if t > threshold_db - 10:
656                         heatmap[f, t] = 1
657                     else:
658                         heatmap[f, t] = 0
659
660         # Create a heatmap of the audiogram
661         heatmap = np.zeros((8000, 100), dtype=np.uint8)
662         for f in range(0, 8000):
663             for t in range(0, 100):
664                 if f < 1000:
665                     if t > threshold_db:
666                         heatmap[f, t] = 1
667                 else:
668                     heatmap[f, t] = 0
669             else:
670                 if f < 2000:
671                     if t > threshold_db + 10:
672                         heatmap[f, t] = 1
673                     else:
674                         heatmap[f, t] = 0
675                 else:
676                     if t > threshold_db - 10:
677                         heatmap[f, t] = 1
678                     else:
679                         heatmap[f, t] = 0
680
681         # Create a heatmap of the audiogram
682         heatmap = np.zeros((8000, 100), dtype=np.uint8)
683         for f in range(0, 8000):
684             for t in range(0, 100):
685                 if f < 1000:
686                     if t > threshold_db:
687                         heatmap[f, t] = 1
688                 else:
689                     heatmap[f, t] = 0
690             else:
691                 if f < 2000:
692                     if t > threshold_db + 10:
693                         heatmap[f, t] = 1
694                     else:
695                         heatmap[f, t] = 0
696                 else:
697                     if t > threshold_db - 10:
698                         heatmap[f, t] = 1
699                     else:
700                         heatmap[f, t] = 0
701
702         # Create a heatmap of the audiogram
703         heatmap = np.zeros((8000, 100), dtype=np.uint8)
704         for f in range(0, 8000):
705             for t in range(0, 100):
706                 if f < 1000:
707                     if t > threshold_db:
708                         heatmap[f, t] = 1
709                 else:
710                     heatmap[f, t] = 0
711             else:
712                 if f < 2000:
713                     if t > threshold_db + 10:
714                         heatmap[f, t] = 1
715                     else:
716                         heatmap[f, t] = 0
717                 else:
718                     if t > threshold_db - 10:
719                         heatmap[f, t] = 1
720                     else:
721                         heatmap[f, t] = 0
722
723         # Create a heatmap of the audiogram
724         heatmap = np.zeros((8000, 100), dtype=np.uint8)
725         for f in range(0, 8000):
726             for t in range(0, 100):
727                 if f < 1000:
728                     if t > threshold_db:
729                         heatmap[f, t] = 1
730                 else:
731                     heatmap[f, t] = 0
732             else:
733                 if f < 2000:
734                     if t > threshold_db + 10:
735                         heatmap[f, t] = 1
736                     else:
737                         heatmap[f, t] = 0
738                 else:
739                     if t > threshold_db - 10:
740                         heatmap[f, t] = 1
741                     else:
742                         heatmap[f, t] = 0
743
744         # Create a heatmap of the audiogram
745         heatmap = np.zeros((8000, 100), dtype=np.uint8)
746         for f in range(0, 8000):
747             for t in range(0, 100):
748                 if f < 1000:
749                     if t > threshold_db:
750                         heatmap[f, t] = 1
751                 else:
752                     heatmap[f, t] = 0
753             else:
754                 if f < 2000:
755                     if t > threshold_db + 10:
756                         heatmap[f, t] = 1
757                     else:
758                         heatmap[f, t] = 0
759                 else:
760                     if t > threshold_db - 10:
761                         heatmap[f, t] = 1
762                     else:
763                         heatmap[f, t] = 0
764
765         # Create a heatmap of the audiogram
766         heatmap = np.zeros((8000, 100), dtype=np.uint8)
767         for f in range(0, 8000):
768             for t in range(0, 100):
769                 if f < 1000:
770                     if t > threshold_db:
771                         heatmap[f, t] = 1
772                 else:
773                     heatmap[f, t] = 0
774             else:
775                 if f < 2000:
776                     if t > threshold_db + 10:
777                         heatmap[f, t] = 1
778                     else:
779                         heatmap[f, t] = 0
780                 else:
781                     if t > threshold_db - 10:
782                         heatmap[f, t] = 1
783                     else:
784                         heatmap[f, t] = 0
785
786         # Create a heatmap of the audiogram
787         heatmap = np.zeros((8000, 100), dtype=np.uint8)
788         for f in range(0, 8000):
789             for t in range(0, 100):
790                 if f < 1000:
791                     if t > threshold_db:
792                         heatmap[f, t] = 1
793                 else:
794                     heatmap[f, t] = 0
795             else:
796                 if f < 2000:
797                     if t > threshold_db + 10:
798                         heatmap[f, t] = 1
799                     else:
800                         heatmap[f, t] = 0
801                 else:
802                     if t > threshold_db - 10:
803                         heatmap[f, t] = 1
804                     else:
805                         heatmap[f, t] = 0
806
807         # Create a heatmap of the audiogram
808         heatmap = np.zeros((8000, 100), dtype=np.uint8)
809         for f in range(0, 8000):
810             for t in range(0, 100):
811                 if f < 1000:
812                     if t > threshold_db:
813                         heatmap[f, t] = 1
814                 else:
815                     heatmap[f, t] = 0
816             else:
817                 if f < 2000:
818                     if t > threshold_db + 10:
819                         heatmap[f, t] = 1
820                     else:
821                         heatmap[f, t] = 0
822                 else:
823                     if t > threshold_db - 10:
824                         heatmap[f, t] = 1
825                     else:
826                         heatmap[f, t] = 0
827
828         # Create a heatmap of the audiogram
829         heatmap = np.zeros((8000, 100), dtype=np.uint8)
830         for f in range(0, 8000):
831             for t in range(0, 100):
832                 if f < 1000:
833                     if t > threshold_db:
834                         heatmap[f, t] = 1
835                 else:
836                     heatmap[f, t] = 0
837             else:
838                 if f < 2000:
839                     if t > threshold_db + 10:
840                         heatmap[f, t] = 1
841                     else:
842                         heatmap[f, t] = 0
843                 else:
844                     if t > threshold_db - 10:
845                         heatmap[f, t] = 1
846                     else:
847                         heatmap[f, t] = 0
848
849         # Create a heatmap of the audiogram
850         heatmap = np.zeros((8000, 100), dtype=np.uint8)
851         for f in range(0, 8000):
852             for t in range(0, 100):
853                 if f < 1000:
854                     if t > threshold_db:
855                         heatmap[f, t] = 1
856                 else:
857                     heatmap[f, t] = 0
858             else:
859                 if f < 2000:
860                     if t > threshold_db + 10:
861                         heatmap[f, t] = 1
862                     else:
863                         heatmap[f, t] = 0
864                 else:
865                     if t > threshold_db - 10:
866                         heatmap[f, t] = 1
867                     else:
868                         heatmap[f, t] = 0
869
870         # Create a heatmap of the audiogram
871         heatmap = np.zeros((8000, 100), dtype=np.uint8)
872         for f in range(0, 8000):
873             for t in range(0, 100):
874                 if f < 1000:
875                     if t > threshold_db:
876                         heatmap[f, t] = 1
877                 else:
878                     heatmap[f, t] = 0
879             else:
880                 if f < 2000:
881                     if t > threshold_db + 1
```

Portable Self Assessment Audiometer - right ear

Sl. No.	Pitch (Frequency Hz)	Hearing Level (Volume dB)	Hearing Loss Range
1	125	0	Normal Hearing (0-15 dB)
2	250	0	Normal Hearing (0-15 dB)
3	500	10	Normal Hearing (0-15 dB)
4	1000	0	Normal Hearing (0-15 dB)
5	2000	0	Normal Hearing (0-15 dB)
6	4000	0	Normal Hearing (0-15 dB)
7	8000	10	Normal Hearing (0-15 dB)

Table 3: Result Analysis



Fig 10.4.6. Portable Self Assessment Audiometer

CHAPTER – 11

CONCLUSION

The project has resulted in a significant achievement: the development of a portable audiometer capable of autonomously and repeatedly diagnosing hearing loss. What distinguishes this portable device is its innovative incorporation of information technology into audiology equipment.

By utilizing the widely-adopted Python programming language, the system is accessible even to individuals without specialized expertise, and Python's user-friendly nature facilitates smooth integration with other devices.

The implementation of automated screening procedures in Python not only cuts down on operational expenses but also enhances the portability of the hearing screening tool. Furthermore, the created audiometer functions as a self-assessment system for hearing, eliminating the requirement for expert operation.

Looking forward, future advancements could focus on integrating Active Noise Cancelling technology into the headphones. This upgrade would effectively reduce ambient noise levels, creating a more conducive environment for conducting hearing tests.

CHAPTER – 12

BIBLIOGRAPHY

- [1] "Ritu Rani; H.T. Patil", "Portable audiometer for detecting hearing disorder at an early stage for cancer patient", "2016 International Conference on Automatic Control and Dynamic Optimization Techniques".
- [2] "Marwa Gargouri; Mondher Chaoui", "Development of hearing self-assessment pure tone audiometer, "2020 IEEE International Conference on Design & Test of Integrated Micro & Nano-Systems (DTS)"
- [3] "M.Dharani kumar Chowdhary, Dr. C. Nagaraja ", "A novel raspberry pi-3 based pure tone audiometer and verification of calibration with standard system, "Journal of Data Acquisition and Processing(2020)".
- [4] "Silvia Figueira; Kevin Nguyen; Shweta Panditrao", "HearThat? - An app for diagnosing hearing loss", "IEEE Global Humanitarian Technology Conference (GHTC 2019)".
- [5] Weinstein BE. Geriatric audiology. New York: Thieme Medical Publisher, Inc. 2000.
- [6] Zeit K. Collaborating with clients & the media to improve public awareness. Adult hearing screening. 2007.
http://www.asha.org/Events/convention/handouts/2007/0413_Zeit_Katrina_3/ Accessed 9 January 2012.
- [7] "Deafness and hearing loss," World Health Organization, Mar. 20, 2019.[Online]. Available: <http://www.who.int/en/news-room/factsheets/detail/deafness-and-hearing-loss>
- [8] American Speech-Language-Hearing Association. Causes of hearing loss in adults 2012. http://www.asha.org/public/hearing/disorders/causes_adults.htm Accessed 3 January
- [9] Demers K. Hearing screening in older adults: a brief hearing loss screener. The Hartford Institute for Geriatric Nursing, New York University College of Nursing 2007.

- [10] Hawthorne G, Hogan A, Giles E, et al. Evaluating the health-related quality of life effects of cochlear implants: A prospective study of an adult cochlear implant program. *Int J Audiol* 2004; 43: 183 – 192.
- [11] R. Filipo et al., “Hyperbaric oxygen therapy with short duration intratympanic steroid therapy for sudden hearing loss,” *Acta OtoLaryngologica*, vol. 132, pp. 475–481, 2012. S. Rajkumar, S. Muttan, and B. Pillai, “Adaptive expert system for audiologists,” in Proc. Int. Conf. Commun. Sig. Process., Feb. 10–22, 2011, pp. 305–309.
- [12] VENCOVSKÝ V. et RUND, F. Pure tone audiometer. 20th Annual Conference Proceeding's Technical Computing. 2012. p. 1- 5.
- [13] FRANKS, John R. Hearing measurement. *Occupational Exposure to Noise: Evaluation, Prevention and Control*. Geneva: World Health Organisation, 2001, p. 183-231
- [14] N. Nakamura, “Development of mobileaudiometer for screening using mobile phones,” in Proc. 26th Annu. Int. Conf. IEEE Eng. Medicine Biol. Soc., Sep. 1–5, 2004, pp. 3369–3372.
- [15] P. G. Jacobs et al., “Development and evaluation of a portable audiometer for high-frequency screening of hearing loss from ototoxicity in homes/clinics,” *IEEE Trans. Biomed. Eng.*, vol. 59, no. 11, pp. 3097–3103, Jul. 2012.
- [16] GAN, Kok Beng, AZEEZ, Dhifaf, UMAT, Cila, et al. Development of a computer-based automated pure tone hearing screening device: a preliminary clinical trial. *Biomedizinische Technik/Biomedical Engineering*, 2012, vol. 57, no 5, p. 323-332.
- [17] RITU RANI, H.T. PATIL “Development and evaluation of a portable audiometer with remote health care”. *International Journal of Industrial Electronics and Electrical Engineering*, ISSN: 2347-698. Volume-4, Issue-6, Jun.-2016
- [18] Mahalakshmi.A, Mohanavalli.M, “PC based audiometer generating audiogram to assess acoustic threshold,” *International Journal of Pure and Applied Mathematics* Vol. 10, no. 12 2018, pp.13939-13944
- [19] Kaplan DM, Shipp DB, Chen JM, Ng AHC, Nedzelski JM. Earlydeafened adult cochlear implant users: assessment of outcomes. *J Otolaryngol* 2003; 32: 245 – 249.

CHAPTER 13

APPENDICES

11.1. SOURCE CODE (PYTHON CODE):

```
import argparse
from datetime import datetime, timedelta
from time import sleep
import numpy as np
import pyaudio
import threading
import pandas as pd
import matplotlib.pyplot as plt
from matplotlib.backends.backend_tkagg import FigureCanvasTkAgg
from pynput.mouse import Listener, Button
import tkinter as tk

# Use interactive backend for displaying plots in a separate window
plt.switch_backend('TkAgg') # You may need to install TkAgg backend if not already installed

class HearingTest:
    def __init__(self):
        self.signal = None
        self.right_data = []
        self.detected = False
        self.start_time = None

    def display_instructions(self):
        """Display instructions in a new window"""
        instructions_window = tk.Tk()
        instructions_window.title("Instructions")
        instructions_label = tk.Label(instructions_window, text="Portable Self Assessment Audiometer", font=("Arial", 16, "bold"))
        instructions_label.pack(padx=10, pady=10)
```

instructions_text = "1. Put on your headphones.\n\n2. Click the left mouse button when you hear pulsing sounds.\n\n3. The test will run for the right ear, playing sounds at different frequencies and volumes."

```
instructions_text_label = tk.Label(instructions_window, text=instructions_text, font=("Arial", 12), justify=tk.LEFT)
instructions_text_label.pack(padx=10, pady=10)
start_button = tk.Button(instructions_window, text="Start Test", command=self.start_test)
start_button.pack(padx=10, pady=10)
instructions_window.mainloop()

def start_test(self):
    """Start the hearing test"""
    self.display_instructions()
    self.run_test()

def player(self, p, repeat=1, ear='right'):
    """Plays sounds with different frequencies and volume levels"""
    volumes = [0, 10, 20, 30, 40, 50, 60, 70, 80, 90] # Adjusted volume levels in dB
    frequencies = [125, 250, 500, 1000, 2000, 4000, 8000] # Adjusted frequencies in Hz
    # Repeat each frequency based on the provided argument
    frequencies = np.repeat(frequencies, repeat)
    stream = p.open(format=pyaudio.paFloat32,
                     channels=1,
                     rate=44100,
                     output=True)
    sleep(0.1)
    for freq in frequencies:
        self.detected = False
        for vol in volumes:
            print(f'Playing frequency: {freq} Hz at volume: {vol} dB for {ear} ear')
            self.signal = [freq, vol, datetime.now()]
            audio_data = (np.sin(2 * np.pi * np.arange(44100 * 0.5) * freq / 44100)).astype(np.float32)
            audio_data = audio_data * 10**((vol / 20))
            stream.write(audio_data.tobytes())
            sleep(2) # Adding 2-second pause after playing each volume level
            if self.detected:
```

```

break

sleep(2) # Adding 2-second pause after playing each frequency
stream.stop_stream()
stream.close()

def on_click(self, x, y, button, pressed):
    """Callback function for mouse clicks"""

    if button == Button.left and pressed:
        if self.signal:
            d = self.signal + [datetime.now()]
            print(f'Recording event: {d}')
            self.right_data.append(d)
            self.detected = True

    def listener(self):
        """Listens to mouse clicks"""

        with Listener(on_click=self.on_click) as listener:
            listener.join()

    def analyse_results(self, data, ear):
        """Stores and visualizes results"""

        now = datetime.now()
        # Load data to DataFrame
        df = pd.DataFrame(data, columns=['frequency', 'volume', 'played', 'heard'])
        df['reaction_time'] = (df['heard'] - df['played']).dt.microseconds // 1000
        # Create audiogram chart
        audiogram_fig = plt.figure()
        ax1 = audiogram_fig.add_subplot(111)
        ax1.plot(df['frequency'], df['volume'], marker='x', linestyle='-', color='black')
        ax1.set(title=f'Audiogram for {ear} ear', ylim=[90, -10], yticks=[90, 80, 70, 60, 50, 40, 30, 20, 10, 0, -10])
        ax1.grid(True)
        ax1.set_ylabel('Hearing Level in decibels (volume in dB)')
        # Add x-axis ticks and labels at the top of the chart
        ax2 = ax1.twiny()
        ax2.set_xlim(ax1.get_xlim())
        ax2.set_xticks(df['frequency'])


```

```

ax2.set_xticklabels(df['frequency'])
ax2.set_xlabel('Pitch (frequency in Hz)')
ax2.xaxis.tick_top()
# Add colored rows for different hearing loss stages
ax1.axhspan(-10, 15, facecolor='green', alpha=0.3, label='Normal Hearing (0-15 dB)')
ax1.axhspan(16, 25, facecolor='yellow', alpha=0.3, label='Slight Hearing Loss (16-25 dB)')
ax1.axhspan(26, 40, facecolor='orange', alpha=0.3, label='Mild Hearing Loss (26-40 dB)')
ax1.axhspan(41, 55, facecolor='red', alpha=0.3, label='Moderate Hearing Loss (41-55 dB)')
ax1.axhspan(56, 70, facecolor='purple', alpha=0.3, label='Moderately Severe Hearing Loss
(56-70 dB)')
ax1.axhspan(71, 90, facecolor='brown', alpha=0.3, label='Severe Hearing Loss (71-90 dB)')
ax1.axhspan(91, 120, facecolor='black', alpha=0.3, label='Profound Hearing Loss (91 dB or
greater)')ax1.legend()
# Save audiogram chart as image
audiogram_fig.savefig(f'./results_{ear}_{now:%Y%m%d%H%M%S}_audiogram.png')
# Display audiogram chart in a new window
audiogram_window = tk.Tk()
audiogram_window.title(f'Audiogram for {ear} ear')
canvas = FigureCanvasTkAgg(audiogram_fig, master=audiogram_window)
canvas.draw()
canvas.get_tk_widget().pack(side=tk.TOP, fill=tk.BOTH, expand=1)
# Display Excel table in a new window
excel_window = tk.Tk()
excel_window.title(f'Portable Self Assessment Audiometer - {ear} ear')
excel_table = tk.Frame(excel_window)
excel_table.grid(row=0, column=0, padx=10, pady=10)
table_label = tk.Label(excel_table, text="Portable Self Assessment Audiometer",
font=("Arial", 16, "bold"))
table_label.grid(row=0, columnspan=5, sticky="w")
headers = ['Sl. No.', 'Pitch (Frequency Hz)', 'Hearing Level (Volume dB)', 'Hearing Loss
Range']
for i, header in enumerate(headers):
    col_label = tk.Label(excel_table, text=header, font=("Arial", 12, "bold"))
    col_label.grid(row=1, column=i, padx=5, pady=5)

```

```

for i, row in df.iterrows():
    sl_no = tk.Label(excel_table, text=i+1, font=("Arial", 12))
    sl_no.grid(row=i + 2, column=0, padx=5, pady=5, sticky="w")
    freq_label = tk.Label(excel_table, text=row['frequency'], font=("Arial", 12))
    freq_label.grid(row=i + 2, column=1, padx=5, pady=5)
    vol_label = tk.Label(excel_table, text=row['volume'], font=("Arial", 12))
    vol_label.grid(row=i + 2, column=2, padx=5, pady=5)
    range_label = tk.Label(excel_table, text=self.get_hearing_loss_range(row['volume']),
                           font=("Arial", 12))
    range_label.grid(row=i + 2, column=3, padx=5, pady=5)
    excel_window.mainloop()

# Create CSV file
df.to_csv(f'./results_{ear}_{now:%Y%m%d%H%M%S}.csv', index=None)

# Create Excel sheet
df_excel = pd.DataFrame(data, columns=['Sl. No.', 'Pitch (Frequency Hz)', 'Hearing Level
(Volume dB)', 'Hearing Loss Range'])
df_excel['Hearing Loss Range'] = df_excel['Hearing Level (Volume
dB)'].apply(self.get_hearing_loss_range)
df_excel.to_excel(f'./results_{ear}_{now:%Y%m%d%H%M%S}.xlsx', index=None)
print("Audiogram chart, CSV file, and Excel sheet created successfully.")

return df

def get_hearing_loss_range(self, volume):
    """Determines the hearing loss range based on volume level"""
    if volume <= 15:
        return 'Normal Hearing (0-15 dB)'
    elif volume <= 25:
        return 'Slight Hearing Loss (16-25 dB)'
    elif volume <= 40:
        return 'Mild Hearing Loss (26-40 dB)'
    elif volume <= 55:
        return 'Moderate Hearing Loss (41-55 dB)'
    elif volume <= 70:
        return 'Moderately Severe Hearing Loss (56-70 dB)'
    elif volume <= 90:

```

```

return 'Severe Hearing Loss (71-90 dB)'

else:
    return 'Profound Hearing Loss (91 dB or greater)'

def run_test(self):
    parser = argparse.ArgumentParser()
    parser.add_argument('-r', '--repeat', help='Number of times each frequency is repeated',
                        type=int, default=1) # Change default value to 1
    args = parser.parse_args()
    self.start_time = datetime.now()
    p = pyaudio.PyAudio()
    # Start listener
    p2 = threading.Thread(target=self.listener, daemon=True)
    p2.start()
    # Run test for the right ear
    print('Testing right ear...')
    self.player(p, repeat=args.repeat, ear='right')
    # Analyse and visualize results for the right ear
    right_df = self.analyse_results(self.right_data, 'right')
    self.right_data = []
    print('Test is finished. Please check visualizations and files.')
    self.display_date_time_duration()

def display_date_time_duration(self):
    now = datetime.now()
    duration = now - self.start_time
    # Display test information in a new window
    info_window = tk.Tk()
    info_window.title("Test Information")
    date_label = tk.Label(info_window, text=f'Date: {self.start_time.strftime("%Y-%m-%d")}', font=("Arial", 12))
    date_label.pack()
    start_time_label = tk.Label(info_window, text=f'Start Time: {self.start_time.strftime("%H:%M:%S")}', font=("Arial", 12))
    start_time_label.pack()
    duration_label = tk.Label(info_window, text=f'Duration: {duration}', font=("Arial", 12))

```

```
duration_label.pack()  
info_window.mainloop()  
if __name__ == '__main__':  
    test = HearingTest()  
    test.start_test()
```

PLAGIARISM REPORT

My project contains 5% plagiarism, and I verified it using the Turnitin plagiarism checker in the MCET library.

Report _Batch 08.pdf

ORIGINALITY REPORT



PRIMARY SOURCES

1	www.coursehero.com Internet Source	3%
2	www.scilit.net Internet Source	1%
3	Marwa Gargouri, Mondher Chaoui, Patrice Wira. "Development of hearing self-assessment pure tone audiometer", 2020 IEEE International Conference on Design & Test of Integrated Micro & Nano-Systems (DTS), 2020 Publication	1%
4	Submitted to University of Canterbury Student Paper	1%

Submission date: 10-May-2024 11:42AM (UTC+0300)

Submission ID: 2372167370

File name: MINI_project.pdf (2.98M)

Word count: 8905

Character count: 52815